



**BD**

Helping all people  
live healthy lives



# Hypodermic Product Catalog

The Basis for Selection of Your  
Medication Delivery Device Needs

# BD PosiFlush™ Prefilled Syringes

**BD™ Prefilled Flush Syringes Are Available In A Wide Variety of Sizes and Solutions To Meet All Your Flush-Protocol Requirements.**

BD Ref. #	Description	NHRIC Number	Box/Case
<b>BD PosiFlush™ Normal Saline Syringes, all with 10 mL diameter</b>			
306544	3 mL	8290-306544	30/Tray 480/Case
306545	5 mL	8290-306545	30/Tray 480/Case
306456	10 mL	8290-306546	30/Tray 480/Case
306549	3 mL with Blunt Plastic Cannula	8290-306549	30/Box 120/Case
306550	5 mL with Blunt Plastic Cannula	8290-306550	30/Box 120/Case
306551	10 mL with Blunt Plastic Cannula	8290-306551	30/Box 120/Case
<b>BD PosiFlush™ SF Saline Syringe</b>			
306553	10mL Sterile Field Syringe	8290-306553	30/Box 120/Case
<b>BD™ Pre-Filled Normal Saline Syringes</b>			
306500	10 mL Fill in 10 mL Syringe	8290-306500	30/Box 120/Case
306502	5 mL Fill in 10 mL Syringe	8290-306502	30/Box 120/Case
306503	3 mL Fill in 10 mL Syringe	8290-306503	30/Box 120/Case
306504	5 mL Fill in 5 mL Syringe	8290-306504	30/Box 120/Case
306505	3 mL Fill in 5 mL Syringe	8290-306505	30/Box 120/Case
306507	3 mL Fill in 3 mL Syringe	8290-306507	30/Box 120/Case
306508	2 mL Fill in 3 mL Syringe	8290-306508	30/Box 120/Case
306518	10 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306518	30/Box 120/Case
306519	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306519	30/Box 120/Case
306522	5 mL Fill in 5 mL Syringe with Blunt Plastic Cannula	8290-306522	30/Box 120/Case
306524	3 mL Fill in 3 mL Syringe with Blunt Plastic Cannula	8290-306524	30/Box 120/Case
<b>BD PosiFlush™ Heparin Lock Flush Syringes, 10 USP units/mL</b>			
306509	6 mL Fill in 10 mL Syringe	8290-306509	30/Box 120/Case
306510	5 mL Fill in 10 mL Syringe	8290-306510	30/Box 120/Case
306511	5 mL Fill in 5 mL Syringe	8290-306511	30/Box 120/Case
306521	3 mL Fill in 10 mL Syringe	8290-306521	30/Box 120/Case
306512	3 mL Fill in 3 mL Syringe	8290-306512	30/Box 120/Case
306525	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306525	30/Box 120/Case
306527	2 mL Fill in 3 mL Syringe	8290-306527	30/Box 120/Case
<b>BD PosiFlush™ Heparin Lock Flush Syringes, 100 USP units/mL</b>			
306513	5 mL Fill in 10 mL Syringe	8290-306513	30/Box 120/Case
306514	3 mL Fill in 10 mL Syringe	8290-306514	30/Box 120/Case
306515	5 mL Fill in 5 mL Syringe	8290-306515	30/Box 120/Case
306516	3 mL Fill in 5 mL Syringe	8290-306516	30/Box 120/Case
306517	3 mL Fill in 3 mL Syringe	8290-306517	30/Box 120/Case
306528	2 mL Fill in 3 mL Syringe	8290-306528	30/Box 120/Case
306531	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	8290-306531	30/Box 120/Case
<b>BD™ SASH Kits: 2 BD Prefilled Normal Saline Syringes, 1 Prefilled Heparin Lock Flush Syringe</b>			
306537	5 mL Fill in 10 mL Syringe – 100 unit/mL	8290-306537	30/Box 240/Case
<b>BD PosiFlush™ Heparin Lock Flush Syringes, 10 USP units/mL</b>			
306413	3 mL	8290-306413	30/Box 480/Case
306414	5 mL	8290-306414	30/Box 480/Case
<b>BD PosiFlush™ Heparin Lock Flush Syringes, 100 USP units/mL</b>			
306423	3 mL	8290-306423	30/Box 480/Case
306424	5 mL	8290-306424	30/Box 480/Case



All products are Preservative-Free and do not contain natural rubber latex

**NOTE:** Bulk non-sterile syringes in bulk packaging are also available. Contact your BD Distributor or your BD Sales Representative for details.



Interlink is a registered trademark of Baxter International, Inc.  
 LifeShield is a registered trademark of Hospina, Inc.  
 SafeLine is a registered trademark of B. Braun Medical, Inc.  
 BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2012 BD

Printed in. USA 07/12 MSS0158

1 Becton Drive  
 Franklin Lakes, NJ 07417  
[www.bd.com](http://www.bd.com)

2771 Bristol Circle  
 Oakville, Ontario L6H 6R5

# BD Eclipse™ Needles and Syringes

## BD Eclipse™ Needles

100 per box, 1,200 per case

305767	25 G x 1 1/2 in.
305761	25 G x 1 in.
305762	23 G x 1 in.
305768	22 G x 1 in.
305763	22 G x 1 1/2 in.
305765	21 G x 1 1/2 in. Thin Wall
305764	21 G x 1 in. Thin Wall
305766	18 G x 1 1/2 in.



## BD Eclipse™ Needles for Luer Lock Syringe

100 per box, 1,200 per case

305759	25 G x 5/8 in.
305758	27 G x 1/2 in.
305757	30 G x 1/2 in.

## BD Eclipse™ Needles for Luer Slip Syringes, featuring SmartSlip™ Technology

100 per box, 1,200 per case

305760	25 G x 5/8 in.
305770	27 G x 1/2 in.
305771	30 G x 1/2 in.

## 1 mL BD Luer Slip Syringe with Detachable Needle

50 per box, 300 per case

305775	30 G x 1/2 in.
305776	27 G x 1/2 in.
305777	25 G x 5/8 in.

## BD Eclipse™ Needles with Syringes

50 per box, 300 per case

### 1 mL BD Luer-Lok™ Syringe with Detachable Needle

305778	30 G x 1/2 in.
305789	27 G x 1/2 in.
305780	25 G x 5/8 in.

### 3 mL BD Luer-Lok™ Syringe with Detachable Needle


305781	25 G x 5/8 in.
305787	25 G x 1 in.
305782	23 G x 1 in.
305788	21 G x 1 1/2 in.
305783	22 G x 1 1/2 in.
305779	21 G x 1 in.
305784	21 G x 1 1/2 in.

### 5 mL BD Luer-Lok™ Syringe with Detachable Needle

305785	22 G x 1 1/2 in.
--------	------------------

### 10 mL BD Luer-Lok™ Syringe with Detachable Needle

305786	22 G x 1 1/2 in.
--------	------------------



# BD Integra™ Needles and Syringes

## BD Integra™ 3 mL Syringe with Detachable Needle

100 per box, 400 per case

305269	3 mL Syringe with 25 G x 5/8 in. Needle
305270	3 mL Syringe with 25 G x 1 in. Needle
305271	3 mL Syringe with 23 G x 1 in. Needle
305272	3 mL Syringe with 22 G x 1 1/2 in. Needle
305273	3 mL Syringe with 21 G x 1 in. Needle
305274	3 mL Syringe with 21 G x 1 1/2 in. Needle



## BD Integra™ 3 mL Syringe Only

100 per box, 400 per case

305283	3 mL Syringe Only
--------	-------------------

## BD Integra™ Needles

100 per box, 1,000 per case

305310	25 G x 5/8 in. Retracting Needle
305311	25 G x 1 in. Retracting Needle
305312	23 G x 1 in. Retracting Needle
305313	22 G x 1 1/2 in. Retracting Needle
305275 <sup>1</sup>	18 G x 1 1/2 in. Blunt Fill Needle
305833 <sup>1,2</sup>	18 G x 1 1/2 in. Blunt Filter Needle

1 – 18 G Blunt Fill and Blunt Filter Needles have a non-retracting needle  
2 – Packaged 50 Box/500 Case



# BD SafetyGlide™ Syringes and Needles

## BD SafetyGlide™ Hypodermic Needle

50 per box, 500 per case

- 305921 27 G x 5/8 in. Needle
- 305901 25 G x 5/8 in. Needle
- 305916 25 G x 1 in. Needle
- 305902 23 G x 1 in. Needle
- 305900 22 G x 1 1/2 in. Needle
- 305915 21 G x 1 in. Needle Thin Wall
- 305917 21 G x 1 1/2 in. Needle Thin Wall
- 305918 18 G x 1 1/2 in. Needle



## BD SafetyGlide™ Insulin Syringe with Permanently Attached Needle

100 per box, 400 per case

- 305930 1 mL Insulin Syringe with 29 G x 1/2 in.
- 305932 1/2 mL Insulin Syringe with 29 G x 1/2 in.
- 305935 3/10 mL Insulin Syringe with 29 G x 1/2 in.
- 305934 1/2 mL Insulin Syringe with 30 G x 5/16 in. Thin Wall
- 305937 3/10 mL Insulin Syringe with 31 G x 5/16 in. Thin Wall



## BD SafetyGlide™ Needle with Syringes

50 per box, 400 per case

### 1 mL BD Syringe with Detachable Needle

- 305903 25 G x 5/8 in.
- 305927 27 G x 5/8 in.

### 3 mL BD Luer-Lok™ Syringe with Detachable Needle

- 305904 25 G x 5/8 in.
- 305905 23 G x 1 in.
- 305906 22 G x 1 1/2 in.
- 305909 21 G x 1 1/2 in. Thin Wall

### 5 mL BD Luer-Lok™ Syringe with Detachable Needle

- 305907 22 G x 1 1/2 in.

### 10 mL BD Luer-Lok™ Syringe with Detachable Needle

- 305908 22 G x 1 1/2 in.



## BD SafetyGlide™ Tuberculin Syringe with Permanently Attached Needle

100 per box, 400 per case

- 305945 1 mL Tuberculin Syringe with 27 G x 1/2 in. Regular Bevel
- 305946 1 mL Tuberculin Syringe with 26 G x 3/8 in. Intradermal Bevel



## BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1000 per case

- 305950 1 mL Allergy Tray with 27 G x 1/2 in. Regular Bevel
- 305951 1 mL Allergy Tray with 26 G x 3/8 in. Intradermal Bevel



# Interlink® System and BD™ Cannulas


## Interlink® Locking Cannula Products for use with Interlink System

100 per box, 1000 per case

303369 Threaded Lock Cannula



303370 Lever Lock Cannula



303365 Y Lock Cannula

## BD Blood Collection Assembly for use with Baxter Interlink®, Hospira LifeShield®, or McGaw SafeLine® Systems

25 per box, 200 per case

303380 BD Vacutainer® Tube Holder and Luer Adapter with Preattached BD Blunt Plastic Cannula



## Syringes with BD™ Cannula for use with Baxter Interlink®, Hospira LifeShield®, or McGaw SafeLine® Systems

100 per box, 800 per case

303391 3 mL Syringe with BD Twinpak™ Cannula

303346 3 mL Syringe with BD Blunt Plastic Cannula

303401 3 mL Syringe with Interlink® Vial Access Cannula

100 per box, 400 per case

303347 5 mL Syringe with BD Blunt Plastic Cannula

303348 10 mL Syringe with BD Blunt Plastic Cannula

303392 5 mL Syringe with BD Twinpak™ Cannula

303393 10 mL Syringe with BD Twinpak™ Cannula

303403 5 mL Syringe with Interlink® Vial Access Cannula

303405 10 mL Syringe with Interlink® Vial Access Cannula



## BD Cannula Products for use with Baxter Interlink®, Hospira LifeShield®, or McGaw SafeLine® Systems

100 per box, 1000 per case

303345 BD Blunt Plastic Cannula

303390 BD Twinpak™ Dual Cannula Device

303367 Interlink® Vial Access Cannula



# BD Safety-Lok™ Safety Syringes

## 1 mL BD Safety-Lok™ Insulin Syringe with Permanently Attached Needle

100 per box, 500 per case

329464 29 G x 1/2 in. U-100 BD Ultra-Fine™

## 1 mL BD Safety-Lok™ Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case

305553 27 G x 1/2 in.

305554 25 G x 5/8 in.



## 3 mL BD Safety-Lok™ Syringe with Detachable Needle

100 per box, 800 per case

309592 25 G x 5/8 in.      309596 22 G x 1 in.

309594 23 G x 1 in.      309595 21 G x 1 1/2 in.

309593 22 G x 1 1/2 in.



## 5 mL BD Safety-Lok™ Syringe with Detachable Needle

50 per box, 400 per case

305561 21 G x 1 1/2 in.



## 10 mL BD Safety-Lok™ Syringe with Detachable Needle

50 per box, 400 per case

305564 21 G x 1 1/2 in.



## BD Safety-Lok™ Syringes only, BD Luer-Lok™ Tip

309606 3 mL Syringe Only      100 box/800 case

305558 5 mL Syringe Only      50 box/400 case

305559 10 mL Syringe Only      50 box/400 case

# BD™ Disposal Solutions



## BD Patient Room Collectors

Available in multiple colors

<b>Recykleen™</b>		
305056	BD Recykleen™ Pearl, 5.4 qt.	20/Case
305057	BD Recykleen™ Red, 5.4 qt.	20/Case
305053	BD Recykleen™ Pearl, 3 gal	10/Case
305054	BD Recykleen™ Red, 3 gal.	10/Case
<b>Horizontal Entry</b>		
305543	Pearl, Nestable, 5.4 qt.	20/Case
305428	Pearl, Nestable, 5.4 qt.	12/Case
305423	Pearl, Nestable, 3 gal.	10/Case
305422	Pearl, Nestable, 2 gal.	10/Case
<b>Side Entry</b>		
305444	Pearl, Nestable, 5.4 qt.	20/Case
305425	Pearl, Nestable, 5.4 qt.	12/Case
305550	Beige, BD Wallmate™, 5.4 qt.	1/Case

## BD Phlebotomy Collectors

305487	1.5 qt.	36/Case
305635	1.0 qt.	60/Case

## BD Multi-Use Nestable Collectors

305480	Clear Top, Large Funnel, 14 qt.	20/Case
305464	Clear Top, 60 cc Funnel, 14 qt.	20/Case
305456	Clear Top, Open, 14 qt.	20/Case
305646	Clear Top, Small Open, 14 qt.	20/Case
305457	Clear Top, Open, 6 gal.	12/Case
305465	Clear Top, Large Funnel, 6 gal.	12/Case

## BD Large and Extra Large Collectors

305024	BD Recykleen™ Hinge Top, 9 gal.	8/Case
305025	BD Recykleen™ Slide Top, 9 gal.	8/Case
305023	BD Recykleen™ Hinge Top, 17 gal	5/Case
305026	BD Recykleen™ Slide Top, 19 gal.	5/Case
305615	Red, Clear Hinged Top, 9 gal.	8/Case
305616	Red, Clear Slide Top, 9 gal.	8/Case
305601	Red, Clear Hinged Top, Gasketed, 9 gal.	8/Case
305602	Red, Clear Slide Top, Gasketed, 9 gal.	8/Case
305610	Red, Clear Hinged Top, 17 gal.	5/Case
305665	Red, Clear Hinged Top, Gasketed, 17 gal.	5/Case
305609	Red, Clear Slide Top, 19 gal.	5/Case
305666	Red, Clear Slide Top, Gasketed, 19 gal.	5/Case

## BD Chemotherapy Collectors

305492	Yellow, Plug Cap, 3 gal.	12/Case
305493	Yellow, Plug Cap, 5 gal.	8/Case
305578	Yellow, Screw Cap, 5 gal.	8/Case
305603	Yellow, Clear Hinged Top, 9 gal.	8/Case
305604	Yellow, Clear Slide Top, 9 gal.	8/Case
305613	Yellow, Clear Slide Top, 19 gal.	5/Case
305614	Yellow, Clear Hinged Top, 17 gal.	5/Case

## BD Multi-Use One-Piece Collectors

305488	Small, 3.3 qt.	24/Case
305489	Medium, 6.9 qt.	12/Case
300467	Medium, Non-Vented, 6.9 qt.	12/Case
305490	Large, 8.2 qt.	12/Case
300470	Large, Non-Vented, 8.2 qt.	12/Case
305491	Open Top, 5 gal.	8/Case
300473	Open Top, Non-Vented, 5 gal.	8/Case
305477	Large Funnel, 5 gal.	8/Case
305577	Screw Cap, 5 gal.	8/Case
305557	Regular Funnel, 1.4 qt.	36/Case

## BD Pharmaceutical Waste Sharps Collectors

305058	BD Recykleen™, 5.4 qt.	20/Case
305055	BD Recykleen™, 3 gal.	10/Case
305330	Counterbalanced door, 5.4 qt.	20/Case
305622	Counterbalanced door, 3 gal.	10/Case
305633	Blue Plug Cap, 3 gal.	12/Case
305634	Blue Hinged Top, 9 gal.	8/Case
305329	Slide Top, 9 gal.	8/Case
305028	Hinge Top, 17 gal.	5/Case
305328	Slide Top, 19 gal.	5/Case

## BD RCRA Hazardous Waste Collectors

<b>DOT</b>		
305073	DOT Port Top, 13 gal.	6/Case
305074	DOT Port Top, 16 gal.	6/Case
<b>Non-DOT</b>		
305066	Plug Cap, 3 gal.	12/Case
305067	Plug Cap, 5 gal.	8/Case
305068	Clear Hinge Top, 9 gal.	8/Case
305069	Clear Slide Top, 9 gal.	8/Case
305071	Clear Hinge Top, 17 gal.	5/Case
305072	Clear Slide Top, 19 gal.	5/Case

## BD Disposal Solutions Cabinets and Accessories

305096	BD Recykleen™ Brown Cabinet, 5.4 qt.	1/Case
305098	BD Recykleen™ Green Cabinet, 5.4 qt.	1/Case
305097	BD Recykleen™ Brown Cabinet, 2 and 3 gal.	1/Case
305099	BD Recykleen™ Green Cabinet, 2 and 3 gal.	1/Case
305017	Locking Wall Cabinet, 5.4 qt.	1/Case
305420	Locking Wall Cabinet, 2 and 3 gal.	1/Case
305475	Nestable, Vertical Entry Wall Cabinet, 3.2 qt.	12/Case
305409	Locking Wall Bracket, 2 and 3 gal.	10/Case
305447	Locking Wall Bracket, 5.4 qt. Collectors	10/Case
305451	One-Piece, Cup Bracket, 1.4 qt.	12/Case
305485	Nestable, Bracket, 8 qt., 14 qt. and 6 gal.	12/Case
305494	Security Lock	12/Case
305496	One-Piece, Bracket, 6.9 qt. and 8.2 qt.	12/Case
305497	One-Piece, Strap, 3.3 qt. through 5 gal.	12/Case
305552	Stabilizer, 2 and 3 gal.	5/Case
305570	Nestable, Stabilizer, 8 qt.	5/Case
305588	Wall/Cart, 2 and 3 gal.	6/Case
305589	Nestable Wire Bracket, 8 qt.	4/Case
305021	Stabilizer	5/Case
305494	Padlock	12/Case
305495	XL Security Bracket	12/Case
305970	Cup Bracket	12/Case

See BD Disposal Solutions Catalog for a complete listing of products available.

# BD PrecisionGlide™ Needles (sterile, single-use)

## BD™ Regular Bevel Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305106	30 G x 1/2 in.	
305128	30 G x 1 in.	
305109	27 G x 1/2 in.	
305136	27 G x 1 1/4 in.	
301629	27 G x 1 1/2 in.	
305110	26 G x 3/8 in.	Intradermal Bevel
305111	26 G x 1/2 in.	
305122	25 G x 5/8 in.	
305124	25 G x 7/8 in.	
305125	25 G x 1 in.	
305127	25 G x 1 1/2 in.	
305143	23 G x 3/4 in.	
305145	23 G x 1 in.	
305120	23 G x 1 1/4 in.	
305155	22 G x 1 in.	
305156	22 G x 1 1/2 in.	
305165	21 G x 1 in.	
305167	21 G x 1 1/2 in.	
305129	21 G x 2 in.	
305175	20 G x 1 in.	
305176	20 G x 1 1/2 in.	
305186	19 G x 1 in.	Thin Wall
305187	19 G x 1 1/2 in.	Thin Wall
305195	18 G x 1 in.	
305196	18 G x 1 1/2 in.	
305197	16 G x 1 in.	
305198	16 G x 1 1/2 in.	
305115	26 G x 5/8 in.	Sub-Q
305193	23 G x 1 in.	Thin Wall IM
305194	23 G x 1 1/2 in.	Thin Wall IM
305185	18 G x 1 1/2 in.	Thin Wall Fill
305177	21 G x 1 in.	Thin Wall IV
305190	21 G x 1 1/2 in.	Thin Wall IV

## BD™ Blunt Fill Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305180	18 G x 1 1/2 in.	
305211	18 G x 1 1/2 in.	Filter (5 Micron)

## BD™ Short Bevel Needles


100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305159	22 G x 1 1/2 in.	
305178	20 G x 1 in.	
305179	20 G x 1 1/2 in.	
305188	19 G x 1 in.	Thin Wall
305189	19 G x 1 1/2 in.	Thin Wall
305199	18 G x 1 1/2 in.	

## Types of needle walls and bevels

The Thin Wall needle is thinner than a Regular Wall needle. The inner diameter is bigger without any increase in the outside diameter. For example, a 21 G Thin Wall needle has the same outside diameter as a 21 G Regular Wall, but has the same inner diameter as a 20 G needle. Thin Wall allows the same flow volume as a Regular Wall needle one gauge larger. Thin Wall is available in 18, 19, 21, 22, and 23 gauge sizes.

### Regular Wall




### Thin Wall

### Regular Bevel

### Short Bevel

### Intradermal Bevel




## Types of syringe tips

### BD Luer-Lok™


#### Tip Sizes:

1 mL  
3 mL  
5 mL  
10 mL  
20 mL  
30 mL  
60 mL



### Slip Tip Sizes:


1 mL  
3 mL  
5 mL  
10 mL  
20 mL  
30 mL  
60 mL



### Eccentric

#### Tip Sizes:

10 mL  
20 mL  
60 mL



### Catheter Tip

#### Sizes:

2 oz. only



BD Conventional  
Needle & Syringe  
Color Chart

30 G

27 G

26 G

25 G

24 G

23 G

22 G

21 G

20 G

19 G

18 G

16 G

# BD™ Syringes & Needles (sterile, single use)

## 10 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 400 per case 1/5 mL

- 309640 22 G x 1 1/2 in.
- 309642 21 G x 1 in.
- 309643 21 G x 1 1/2 in.
- 309644 20 G x 1 in.
- 309645 20 G x 1 1/2 in.

100 per box, 800 per case

- 305064 18 G x 1 1/2 in. BD™ Blunt Fill Needle

### BD Luer-Lok™ Tip, Thin Wall

100 per box, 400 per case 1/5 mL

- 309639 21 G x 1 1/2 in.



## 10 mL BD Syringe Only

100 per box, 400 per case 1/5 mL

- 309604 BD Luer-Lok™ Tip
- 301604 Slip Tip
- 305482 Eccentric Tip

25 per box, 100 per case

- 309695 BD Luer-Lok™ Tip Control Syringe

100 per box, 400 per case

- 305462 Lab Syringe Eccentric Tip Non-Sterile, Clean, Ready-to-use Bulk-Pak



## 20 mL BD Syringe Only

48 per box, 192 per case 1 mL, 3/4 oz. in 1/8 oz.

- 302830 BD Luer-Lok™ Tip
- 302831 Slip Tip

120 per box, 480 per case


- 300613 Eccentric Tip



## 30 mL BD Syringe Only

56 per box, 224 per case 1 mL, 1 oz. in 1/4 oz.

- 302832 BD Luer-Lok™ Tip
- 302833 Slip Tip




## 60 mL BD Syringe Only

40 per box, 160 per case 1 mL, 2 oz. in 1/2 oz.

- 309654 Slip Tip
- 309653 BD Luer-Lok™ Tip

60 per box, 240 per case


- 300866 Eccentric Tip



## 2 OZ BD Catheter Tip Syringe Only (with Tip Shield)

40 per box, 160 per case 1 mL, 2 oz. in 1/2 oz.

- 309620 Catheter Tip



## BD Pharmacy Convenience Paks

### Tray Contents:

- 309701 1 mL Luer Slip Tip 25 syringes only per sterile tray, 12 trays per case of 300
- 309702 3 mL BD Luer-Lok™ Tip 25 syringes only per sterile tray, 12 trays per case of 300
- 309703 5 mL BD Luer-Lok™ Tip 25 syringes only per sterile tray, 12 trays per case of 300
- 309605 10 mLBD Luer-Lok™ Tip 20 syringes only per sterile tray, 12 trays per case of 240
- 305617 20 mLBD Luer-Lok™ Tip 10 each Syringes only per sterile tray, 12 trays per case of 120
- 305618 30 mLBD Luer-Lok™ Tip 10 each Syringes only per sterile tray, 12 trays per case of 120
- 309680 60 mLBD Luer-Lok™ Tip 20 each Syringes only per sterile tray, 6 trays per case of 120



## BD Luer-Lok™ Sterile Syringe Tip Cap

200 per box, 2000 per case

- 305819 BD Syringe Tip Cap



## BD Luer Slip™ Sterile Syringe Tip Cap Tray

50 per tray, 1,000 per box, 4,000 per case

- 305822 Syringe Tip Cap Tray






# BD™ Tuberculin Syringes

## 1/2 mL BD Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case 1/100 mL

305620 27 G x 1/2 in. Regular Bevel



## 1 mL BD Tuberculin Syringe with Detachable Needle, Slip Tip

100 per box, 800 per case 1/100 mL

309623 27 G x 1/2 in.  
309625 26 G x 3/8 in. Intradermal Bevel  
309626 25 G x 5/8 in.  
309624 21 G x 1 in.



## 1 mL BD Tuberculin Syringe Only

200 per box, 1600 per case 1/100 mL

309659 Slip Tip



# BD™ Syringes & Needles (sterile, single-use)

## 1 mL BD Syringe/Needle Combination, Slip Tip

100 per box, 800 per case 1/100 mL

309597 26 G x 5/8 in. Sub-Q



## 1 mL BD Luer-Lok™ Syringe/Needle Combination

100 per box, 800 per case 1/100 mL

309637 20 G x 1 in.



## 1 mL BD Luer-Lok™ Syringe Only

100 per box, 800 per case 1/100 mL

309628 BD Luer-Lok™



## 3 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip™, Thin Wall

100 per box, 800 per case 1/10 mL

309587 26 G x 5/8 in. Sub-Q  
309588 23 G x 1 in. IM  
309589 23 G x 1 1/2 in. IM  
309598 21 G x 1 in. IV  
309599 21 G x 1 1/2 in. IV



## 3 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 800 per case 1/10 mL

309582 25 G x 1 1/2 in.	309577 21 G x 1 1/2 in.
309581 25 G x 1 in.	309575 21 G x 1 in.
309570 25 G x 5/8 in.	309579 20 G x 1 1/2 in.
309571 23 G x 1 in.	309578 20 G x 1 in.
309574 22 G x 1 1/2 in.	309580 18 G x 1 1/2 in.
309573 22 G x 1 1/4 in.	305060 18 G x 1 1/2 in. Blunt Fill Needle
309572 22 G x 1 in.	
309569 22 G x 3/4 in.	



## 3 mL BD Syringe Only

200 per box, 800 per case 1/10 mL

309657 BD Luer-Lok™ Tip  
309656 Slip Tip



## 5 mL BD Syringe/Needle Combination BD Luer-Lok™ Tip

100 per box, 400 per case 1/5 mL

309630 22 G x 1 in.	309634 20 G x 1 in.
309631 22 G x 1 1/2 in.	309635 20 G x 1 1/2 in.
309632 21 G x 1 in.	305062 18 G x 1 1/2 in. Blunt Fill Needle
309633 21 G x 1 1/2 in.	



## 5 mL BD Syringe Only

125 per box, 500 per case 1/5 mL

309646 BD Luer-Lok™ Tip  
309647 Slip Tip



# BD™ Insulin Syringes

## 3/10 mL BD Lo-Dose™ Insulin Syringe with Permanently Attached Needle

Self-Contained Package. 100 per box, 500 per case

328438 31 G x 5/16 in. U-100 BD Ultra-Fine™ Short Needle  
200 per box, 1,600 per case

309301 29 G x 1/2 in. U-100 BD Ultra-Fine



## 1/2 mL BD Lo-Dose Insulin Syringe with Permanently Attached Needle

Blister Package. 100 per box, 500 per case

329461 28 G x 1/2 in. U-100 BD Micro-Fine IV (Orange)

Self-Contained Package. 100 per box, 500 per case

329465 28 G x 1/2 in. U-100 BD Micro-Fine IV (Orange)

329466 29 G x 1/2 in. U-100 BD Ultra-Fine

328468 31 G x 5/16 in. U-100 BD Ultra-Fine Short Needle

200 per box, 1,600 per case

309306 29 G x 1/2 in. U-100 BD Ultra-Fine



NOTE: All of the syringes in this section are in an Easy-Peel Blister Package unless otherwise noted as Self-Contained Package.

## 1 mL BD Insulin Syringe with Permanently Attached Needle

Self-Contained Package. 100 per box, 500 per case

329410 28 G x 1/2 in. U-100 BD Micro-Fine IV (Orange)

329412 27 G x 5/8 in. U-100 BD Micro-Fine IV (Orange)

329411 29 G x 1/2 in. U-100 BD Ultra-Fine

328418 31 G x 5/16 in. U-100 BD Ultra-Fine Short Needle

Self-Contained Package. 200 per box, 1600 per case

309311 29 G x 1/2 in. U-100 BD Ultra-Fine

Blister Package. 100 per box, 500 per case

329424 28 G x 1/2 in. U-100 BD Micro-Fine IV (Orange)

329420 28 G x 1/2 in. U-100 BD Micro-Fine IV  
Single Unit Graduation



## 1 mL BD Insulin Syringe with Detachable Needle

Blister Package. 100 per box, 800 per case

329651 25 G x 5/8 in. U-100 (Orange/Blue)

329652 26 G x 1/2 in. U-100 (Orange/Tan)


329622 25 G x 1 in. U-100 (Orange/Blue)



## 1 mL BD Insulin Syringe Only

Blister Package. 200 per box, 1600 per case

329654 U-100 Slip Tip (Orange)



100 per box, 800 per case

309629 BD Luer-Lok™, 1 mL

# BD™ Allergy Syringes

## BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1,000 per case

305950 1 mL Allergy Tray with  
27 G x 1/2 in. Regular Bevel

305951 1 mL Allergy Tray with  
26 G x 3/8 in. Intradermal Bevel



## 1/2 mL BD Allergist Tray with Permanently Attached Needle

25 syringes per tray, 1,000 per case

305535 27 G x 1/2 in. Regular Bevel

305536 27 G x 3/8 in. Intradermal Bevel



## 1 mL BD Allergist Tray with Permanently Attached Needle

25 syringes per tray, 1,000 per case

305537 26 G x 1/2 in. Regular Bevel

305538 26 G x 1/2 in. Intradermal Bevel

305539 26 G x 3/8 in. Intradermal Bevel

305540 27 G x 1/2 in. Regular Bevel

305541 27 G x 3/8 in. Intradermal Bevel

305542 27 G x 3/8 in. Regular Bevel



## 1 mL BD Allergy Syringe with Permanently Attached Needle

100 per box, 500 per case

305500 28 G x 1/2 in. Regular Bevel

