


ETHICON

**SURGICAL
STAINLESS STEEL
WIRE**

Monofilament and Multifilament



SURGICAL STAINLESS STEEL WIRE

Monofilament and Multifilament

Colour

Silver.

Construction

Multifilament and monofilament available.

Non-Absorbable

Provides permanent tensile strength retention in tissue.

Multifilament Range

Gauge sizes 4-0 and 3-0 (U.S.P.). Supplied as needled sutures.

Monofilament Range








Gauge sizes 1 to 7. Supplied as needled sutures.

Sterilisation


By Irradiation.

These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).




Stainless Steel Wire

Needle length and description		Length		Gauge					
☉ 1/2 circle Taper Point				2-0 3	0 3.5	1 4	2 5	4 6	5 7
 MH-1	31mm	75cm	multi	W311					
▼ 1/2 circle TAPERCUT*				0 3.5	1 4	2 5	4 6	5 7	7 9
 V-37 heavy	40mm	4x45cm	mono		M660G				
STRAIGHT PACK									
 V-40 heavy	48mm	4x45cm	mono				M651G	M650G	
STRAIGHT PACK									
▼ 1/2 circle Reverse Cutting				0 3.5	1 4	2 5	4 6	5 7	7 9
 CPX	48mm	4x45cm	mono						M624G
STRAIGHT PACK									
▼ 1/2 circle Reverse Cutting Rotating Needle				0 3.5	1 4	2 5	4 6	5 7	7 9
 CPXX	55mm	4x45cm	mono					M635G	
STRAIGHT PACK									
 CPXX	55mm	45cm 6 x 45cm 75cm 2 x 75cm	mono mono mono mono					EH7216G MEH7216G W905 W902	
STRAIGHT PACK									
 CPXX	55mm	75cm	mono					W995	

Stainless Steel Wire

Needle length and description	Length			Gauge					
▲ 1/2 circle TROCAR POINT			$\frac{2-0}{3}$	$\frac{0}{3.5}$	$\frac{1}{4}$	$\frac{2}{5}$	$\frac{4}{6}$	$\frac{5}{7}$	$\frac{7}{9}$
TR-55 	55mm	75cm	mono					W945	

Stainless Steel Wire

Needle length and description	Length					
		5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
PULL-OUT TENDON SUTURES						
50mm straight Taper Point  22mm 1/2 circle reverse cutting	30cm Multi with barb				FSN13	
TEMPORARY CARDIAC PACING WIRE <i>insulated with blue polyethylene.</i>						
Incorporating a 60mm straight needle which has a 'breakaway' feature enabling the remaining butt to be used as a connecting plug to the Pacing Unit.						
			Diameter			
			0.3mm			
17mm TAPERCUT* 1/2 circle  60mm straight cutting	60cm Multi		FEP13E			
26mm TAPERCUT* 1/2 circle  60mm straight cutting	60cm Multi		FEP15E			

