



Xenon Light Projector
AUTO LP 5132 **HIGHT**
300 Watt Xenon


Automatic high power light source with
dialog-function

Xenon Light Projector LP 5132

Innovations for Intelligent Illumination

□ Intelligent


- in *dialog* mode* the LP HighLight communicates through a data cable with your endocamera and is controlled optimally according to the OR situation



□ Cost-saving

- by optimum control in *dialog* mode the light energy is always adjusted to the suitable minimum power. This spares your endoscopic equipment and reduces costs by prolonged life
- the new 300 Watt lamp module is more favorable in price, so that follow-on costs for the lamp as consumable article are minimised

□ Communicative

- simple operation by touch screen and voice control in our integrated **CORE** OR solution 

□ Safe

- no glare caused by light output if no light cable has been plugged in

□ Simplified repair

- easy replacement of the xenon lamp module without tools via the side access panel

HighLight light projector set consisting of:

Light projector (5132.002), 300 Watt Xenon Lamp module (2431.121) and cable accessories, power connection
100 - 240 V, 50 / 60 Hz5132.0021

* Only in connection with endocamera with *dialog* capability

Examples of combinations with maximum functional equipment

1 CCD



1 CCD
Camera
system
5520

LP 5132
HighLight

3 CCD



3 CCD
Camera
system
5508
Xtended

LP 5132
HighLight

RICHARD WOLF GmbH · 75434 Knittlingen · PF 1164 · Phone +49 7043 - 35-0 · Fax +49 7043 - 3 53 00 · GERMANY · e-mail: info@richard-wolf.com · www.richard-wolf.com

AUSTRIA:

RICHARD WOLF Austria
Ges.m.b.H.
Wilhelminenstraße 93 a
A-1160 Wien
Tel.: 01-405 51 51
Fax: 01-405 51 51 45
info@richard-wolf.at

BELGIUM:

N.V. Endoscopie
RICHARD WOLF Belgium S.A.
Industriezone Drogen, Landegemstraat 6
B-9031 Gent-Drogen
Tel.: 09.280.81.00
Fax: 09.282.92.16
endoscopy@richard-wolf.be

FRANCE:

RICHARD WOLF France S.A.R.L.
Rue Daniel Berger
Z.A.C. La Neuville
F-51100 Reims
Tel.: 03.26.87.02.89
Fax: 03.26.87.60.33
endoscopes@richardwolf.fr

UK:

RICHARD WOLF UK Ltd
Waterside Way
Wimbledon
SW17 0HB
Tel.: 020 8944 7447
Fax: 020 8944 1311
admin@richardwolf.uk.com

USA:

RICHARD WOLF
Medical Instruments Corp.
353 Corporate Woods Parkway
Vernon Hills, Illinois 60061
Tel.: 847-913 1113
Fax: 847-913 1488
sales&marketing@richardwolfusa.com