

## HEMOCLIP PLUS EASY LOAD

### SIMPLER APPLICATION.


Since Weck pioneered the first clip in 1963 the Hemoclip system has become the worldwide standard.

Today this clip is available in the form of a unique EASY-LOAD clip bank, which stands for increased speed, comfort and dexterity during surgery.

When you add up all the advantages, including the new clip bank, the heart-shaped wire, the complete closure of the clip and the profiled instruments, you arrive at a premium clip system for worldwide use in the operating theatre.

### OUR ADVANTAGES:

- new clip bank
- heart-shaped wire
- complete closure
- profiled clip appliers
- easy loading
- colour coding



EASY LOAD - clip loading made simpler

## HEMOCLIP PLUS TITAN LIGATION CLIPS

Ref.	Size	Tape (yes/no)	Clips/Box	Bank /Box	Clips/Bank
533 835	small	yes	180	18	10
533 837	small	no	180	18	10
534 835	small strongpoint	yes	180	18	10
533 735	small	yes	300	12	25
533 737	small	no	300	12	25
534 735	small strongpoint	yes	300	12	25
533 800	medium	yes	180	18	10
533 802	medium	no	180	18	10
533 700	medium	yes	250	10	25
533 702	medium	no	250	10	25
533 860	medium/large	yes	120	12	10
533 862	medium/large	no	120	12	10
533 870	large	yes	120	12	10
533 872	large	no	120	12	10

## HEMOCLIP PLUS LIGATION CLIP APPLIER

Ref.	Size	Length (cm)	Version
533 140	small	15	curved
533 150	small	20	curved
533 151	small	28	curved
533 152	small	20	right-angled
534 140	small strongpoint	15	curved
534 150	small strongpoint	20	curved
534 151	small strongpoint	28	curved
533 105	medium	20	right-angled
533 106	medium	27	right-angled
533 108	medium	28	45°
533 109	medium	15	curved
533 110	medium	20	curved
533 111	medium	27	curved
533 164	medium/large	28	45°
533 165	medium/large	20	curved
533 166	medium/large	27	curved
533 167	medium/large	20	right-angled
533 168	medium/large	27	right-angled
533 175	large	20	curved
533 177	large	20	right-angled
533 178	large	27	right-angled
533 180	large	27	curved
533 181	large	28	45°



Available  
with or without  
adhesive back

Colour coding:

- small
- small strongpoint
- medium
- medium / large
- large

## HEMOCLIP PLUS LIGATION CLIP SYSTEM ACCESSORIES

Ref.	Designation	Length (cm)	Version
523 120	Removal tongs	16	straight
523 121	Removal tongs	27	angled
523 130	Cartridge	18	

Manufacturer: Weck Closure Systems, Research Triangle Park, NC 27709

Distributed by:

TELEFLEX MEDICAL Europe • IDA Business Park • Athlone • Co Westmeath • Ireland  
Tel. +353 (0)9 06 46 08 00 • Fax +353 (0)14 37 07 73 • orders.intl@teleflexmedical.com

www.teleflexmedical.com