

# IMPLANT CATALOGUE

EFFECTIVE JANUARY 2012



**MENTOR**<sup>®</sup>  
BREAST IMPLANTS



Make life more beautiful™

NOT FOR DISTRIBUTION IN USA  
APPROVED FOR EMEA REGION

© Mentor Worldwide LLC 2012 0907007 Rev C

In 2002, we completed a new plant in Leiden, The Netherlands. At 6,000 square meters, it is the largest of its kind in the world. This plant contains an automated breast implant dipping machine, which produces every implant shell in a precisely uniform manner. To meet an ever-increasing demand for our products, the new plant produces more than 300,000 breast implants a year, distributed to more than 100 nations. At Mentor all our products are made under strict standards of design and testing.



Founded in 1969, Mentor Worldwide LLC is a leading manufacturer and supplier of medical products for the global Aesthetic medicine market and Reconstructive surgery. As the world's leading maker of high-quality breast implants for over 25 years, our experience results in quality products that you can rely on.


Headquartered in Santa Barbara, California, Mentor has manufacturing and research operations in the United States, Mauritius, and the Netherlands. In 1995, Mentor became the first breast implant manufacturer to achieve CE-marking for its products.

Mentor has worked closely with surgeons throughout the world to develop innovative products such as SPECTRA™ Adjustable Gel Breast Implant, BECKER™ Expander/Breast Implants, and the CPG™ Gel Breast Implants.

Mentor's manufacturing capability is exceptional. We have substantial experience in biomaterials and biopolymers manufacturing. At the same time, our regulatory expertise has enabled the Company to bring products to market in the world's most stringently regulated health care environments.

Mentor is a manufacturer with a strong customer focus. Our sales representatives are skilled in technical product issues, enabling them to act as a valuable information resource. For questions or additional information, please contact your local Mentor branch or Johnson & Johnson office (listed on back cover) or distributor.

On Jan 23, 2009, Mentor was acquired by Johnson & Johnson and is a stand-alone business unit reporting through ETHICON Inc., a leading provider of suture, mesh, hemostats and other products for a wide range of surgical procedures. By joining forces, we aspire to be the trusted global leader in aesthetic medicine among both consumers and professionals by providing a broad range of innovative, science and clinical-based solutions to maintain, enhance, and restore self-esteem and quality of life. Simply said, our goal is to make life more beautiful.



## CONTENTS

<b>MemoryGel™ Breast Implants, A More Evolved Implant</b> .....	<b>1</b>
About MemoryGel™ Implants	
MemoryGel™ Implants Product Scale	
<b>All Cohesive, All the Time™</b> .....	<b>2</b>
<b>Round Gel Breast Implants, Cohesive I™</b> .....	<b>3</b>
SILTEX™ Round Breast Implants, Moderate Profile, Cohesive I™	
SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive I™	
SILTEX™ Round Breast Implants, High Profile, Cohesive I™	
SILTEX™ Round Breast Implants, Ultra High Profile, Cohesive I™	
Smooth Round Breast Implants, Low Profile, Cohesive I™	
Smooth Round Breast Implants, Moderate Profile, Cohesive I™	
Smooth Round Breast Implants, Moderate Plus Profile, Cohesive I™	
Smooth Round Breast Implants, High Profile, Cohesive I™	
Smooth Round Breast Implants, Ultra High Profile, Cohesive I™	
<b>Round Gel Breast Implants, Cohesive II™</b> .....	<b>7</b>
SILTEX™ Round Breast Implants, Moderate Profile, Cohesive II™	
SILTEX™ Round Breast Implants, Moderate Plus Profile, Cohesive II™	
SILTEX™ Round Breast Implants, High Profile, Cohesive II™	
<b>CPG™ Gel Breast Implants, Cohesive III™</b> .....	<b>9</b>
CPG™ Implants 331, Cohesive III™, Tall Height, Moderate Projection	
CPG™ Implants 321, Cohesive III™, Medium Height, Moderate Projection	
CPG™ Implants 311, Cohesive III™, Low Height, Moderate Projection	
CPG™ Implants 332, Cohesive III™, Tall Height, Moderate Plus Projection	
CPG™ Implants 322, Cohesive III™, Medium Height, Moderate Plus Projection	
CPG™ Implants 312, Cohesive III™, Low Height, Moderate Plus Projection	
CPG™ Implants 333, Cohesive III™, Tall Height, High Projection	
CPG™ Implants 323, Cohesive III™, Medium Height, High Projection	
CPG™ Implants 313, Cohesive III™, Low Height, High Projection	
CPX2™ Tissue Expanders / CPG™ Breast Implants	
Sizing Guide – Two-Stage Breast Reconstruction	
<b>BECKER™ Expander/Breast Implants</b> .....	<b>15</b>
SILTEX™ Round BECKER™ 25 Expander/Breast Implants, Cohesive I™	
SILTEX™ Round BECKER™ 50 Expander/Breast Implants, Cohesive I™	
SILTEX™ CONTOUR PROFILE™ BECKER™ 35 Expander/Breast Implants, Cohesive II™	
<b>SPECTRA™ Adjustable Gel Breast Implants</b> .....	<b>16</b>
SILTEX™ Round SPECTRA™ Adjustable Gel Implants, Cohesive I™	
Smooth Round SPECTRA™ Adjustable Gel Implants, Cohesive I™	
<b>Sterile Resterilizable Gel Breast Implant Sizers</b> .....	<b>17</b>
Round Breast Implant Sizers, Low Profile	
Round Breast Implant Sizers, Moderate Profile	
Round Breast Implant Sizers, Moderate Plus Profile	
Round Breast Implant Sizers, High Profile	
Round Breast Implant Sizers, Ultra High Profile	

## CPG™ Gel Breast Implant Sizers

CPG™ Implant Sizers 331, Cohesive III™, Moderate Projection	
CPG™ Implant Sizers 321, Cohesive III™, Medium Height, Moderate Projection	
CPG™ Implant Sizers 311, Cohesive III™, Low Height, Moderate Projection	
CPG™ Implant Sizers 332, Cohesive III™, Tall Height, Moderate Plus Projection	
CPG™ Implant Sizers 322, Cohesive III™, Medium Height, Moderate Plus Projection	
CPG™ Implant Sizers 312, Cohesive III™, Low Height, Moderate Plus Projection	
CPG™ Implant Sizers 333, Cohesive III™, Tall Height, High Projection	
CPG™ Implant Sizers 323, Cohesive III™, Medium Height, High Projection	
CPG™ Implant Sizers 313, Cohesive III™, Low Height, High Projection	

## CONTOUR PROFILE™ Tissue Expanders .....

SILTEX™ CONTOUR PROFILE™ Expander, Low Height, Style 6100	<b>23</b>
SILTEX™ CONTOUR PROFILE™ Expander, Medium Height, Style 6200	
SILTEX™ CONTOUR PROFILE™ Expander, Tall Height, Style 6300	

## Other Tissue Expanders .....

Smooth Round Tissue Expanders	<b>24</b>
Smooth Elliptical Tissue Expanders	
Smooth Crescent Tissue Expanders	
Smooth Rectangle Tissue Expanders	

## Round Saline Breast Implants .....

SILTEX™ Round Breast Implants, Moderate Profile	<b>25</b>
Smooth Round Breast Implants, Moderate Profile	
Smooth Round Breast Implants, Moderate Plus Profile	
Smooth Round Breast Implants, High Profile	

## CONTOUR PROFILE™ Saline Breast Implants .....

SILTEX™ CONTOUR PROFILE™ Implant, Moderate Profile	<b>27</b>
SILTEX™ CONTOUR PROFILE™ Implant, High Profile	

## SPECTRUM™ Post-Operatively Adjustable Saline Breast Implants .....

Smooth Round SPECTRUM™ Implants	<b>28</b>
SILTEX™ Round SPECTRUM™ Implants	
SILTEX™ CONTOUR PROFILE™ SPECTRUM™ Implants	

## Sterile Saline Breast Implant Sizers .....

Round Breast Implant Sizers, Moderate Profile	<b>29</b>
Round Breast Implant Sizers, Moderate Plus Profile	
Round Breast Implant Sizers, High Profile	
CONTOUR PROFILE™ Saline Breast Implant Sizers, High Profile	

## Accessories .....

Aseptic Transfer System	<b>31</b>
Injection Port Detector	
Injection Dome	

## Important Safety Information .....

## Contact Information .....

**33**

# MemoryGel™

BREAST IMPLANTS

## A More Evolved Implant

- Strict US and EU manufacturing standards
- Highly compliant shell for ease of placement
- Lifetime product replacement policy
- PatientSafe coverage limited warranty or PatientCare Promise limited warranty†
- Low, Moderate, Moderate Plus, High, and Ultra High projections
- Proprietary silicone gel formula for a natural touch that resembles breast tissue

## About MemoryGel™ Breast Implants


MENTOR® MemoryGel™ Breast Implants feature a proprietary cohesive silicone gel formulation used to fill all MENTOR® Silicone Gel-Filled Breast Implants around the world. By varying the components of the gel, Mentor is able to produce a wide selection of products ranging from very soft to very firm.

### MemoryGel™ Breast Implants Product Scale

Product Profile	Cohesive I™	Cohesive II™	Cohesive III™
Round Low	X		
Round Moderate Plus	X	X	
Round Moderate	X	X	
Round High	X	X	
Round Ultra High	X		
CPG™ Gel Breast Implant			X
Round BECKER™ 25 Expander/Breast Implant	X		
Round BECKER™ 50 Expander/Breast Implant	X		
CONTOUR PROFILE™ BECKER™ 35 Expander/Breast Implant		X	
SPECTRA™ Adjustable Gel Breast Implant	X		

† In case of eligibility, depending on the date of surgery, the implant will fall under one of the above warranties.

## All Cohesive, All the Time.™



Advanced MemoryGel™ Breast Implants.  
The more evolved implant.

All MENTOR® Silicone Gel-Filled Breast Implants contain gel that is cohesive, safe and aesthetically pleasing.

### Cohesive I™

The standard cohesive level gel used in MENTOR® Breast Implants. This is the softest gel and has been preferred for years in textured and smooth Round Low, Moderate, Moderate Plus, High, and Ultra High Profile Gel Breast Implants as well as the Round SPECTRA™ Adjustable Gel Breast Implants and BECKER™ Expander/Breast Implants.

### Cohesive II™

A slightly firmer gel, for surgeons wanting a firmer feeling implant. This gel is used in textured Round Moderate, Moderate Plus, High Profile Gel Breast Implants and CONTOUR PROFILE™ BECKER™ Expander/Breast Implants.

### Cohesive III™

Mentor's most cohesive gel, providing shape retention with a pleasing level of firmness for optimal aesthetic results. This gel is used in the MENTOR® CPG™ Gel Breast Implants.

Round Gel Breast Implants,  
Cohesive I™

Round Gel Breast Implants,  
Cohesive I™ CONT'D

SILTEX™ Round Breast Implants,  
Moderate Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8,8 cm	2,5 cm	354-1007	RSZ-7100
125 cc	9,3 cm	2,8 cm	354-1257	RSZ-7125
150 cc	10,2 cm	2,7 cm	354-1507	RSZ-7150
175 cc	10,7 cm	2,8 cm	354-1757	RSZ-7175
200 cc	11,2 cm	2,8 cm	354-2007	RSZ-7200
225 cc	11,4 cm	3,0 cm	354-2257	RSZ-7225
250 cc	11,5 cm	3,3 cm	354-2507	RSZ-7250
275 cc	12,4 cm	3,4 cm	354-2757	RSZ-7275
300 cc	12,6 cm	3,5 cm	354-3007	RSZ-7300
325 cc	12,9 cm	3,6 cm	354-3257	RSZ-7325
350 cc	13,4 cm	3,7 cm	354-3507	RSZ-7350
375 cc	13,4 cm	3,8 cm	354-3757	RSZ-7375
400 cc	13,5 cm	3,9 cm	354-4007	RSZ-7400
450 cc	13,9 cm	4,1 cm	354-4507	RSZ-7450
500 cc	14,2 cm	4,3 cm	354-5007	RSZ-7500
550 cc	14,8 cm	4,4 cm	354-5507	RSZ-7550
600 cc	15,4 cm	4,5 cm	354-6007	RSZ-7600
700 cc	16,8 cm	4,3 cm	354-7007	RSZ-7700
800 cc	17,2 cm	4,6 cm	354-8007	RSZ-7800

SILTEX™ Round Breast Implants,  
High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
125 cc	8,4 cm	3,6 cm	354-4125	RSZ-1254
150 cc	8,9 cm	3,8 cm	354-4150	RSZ-1504
175 cc	9,4 cm	4,0 cm	354-4175	RSZ-1754
200 cc	9,9 cm	4,1 cm	354-4200	RSZ-2004
225 cc	10,2 cm	4,3 cm	354-4225	RSZ-2254
250 cc	10,5 cm	4,5 cm	354-4250	RSZ-2504
275 cc	10,9 cm	4,6 cm	354-4275	RSZ-2754
300 cc	11,1 cm	4,7 cm	354-4300	RSZ-3004
325 cc	11,5 cm	4,8 cm	354-4325	RSZ-3254
350 cc	11,7 cm	4,9 cm	354-4350	RSZ-3504
375 cc	12,0 cm	5,0 cm	354-4375	RSZ-3754
400 cc	12,3 cm	5,1 cm	354-4400	RSZ-4004
425 cc	12,5 cm	5,2 cm	354-4425	RSZ-4254
450 cc	12,7 cm	5,2 cm	354-4450	RSZ-4504
500 cc	13,2 cm	5,4 cm	354-4500	RSZ-5004
550 cc	13,5 cm	5,6 cm	354-4550	RSZ-5504
600 cc	14,0 cm	5,7 cm	354-4600	RSZ-6004
650 cc	14,3 cm	5,8 cm	354-4650	RSZ-6504
700 cc	14,7 cm	6,0 cm	354-4700	RSZ-7004
800 cc	15,4 cm	6,3 cm	354-4800	RSZ-8004

SILTEX™ Round Breast Implants,  
Moderate Plus Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
100 cc	8,1 cm	2,7 cm	354-1001	RSZ-1001
125 cc	8,8 cm	2,9 cm	354-1251	RSZ-1251
150 cc	9,4 cm	3,0 cm	354-1501	RSZ-1501
175 cc	10,0 cm	3,2 cm	354-1751	RSZ-1751
200 cc	10,5 cm	3,3 cm	354-2001	RSZ-2001
225 cc	10,9 cm	3,5 cm	354-2251	RSZ-2251
250 cc	11,3 cm	3,6 cm	354-2501	RSZ-2501
275 cc	11,7 cm	3,7 cm	354-2751	RSZ-2751
300 cc	12,0 cm	3,7 cm	354-3001	RSZ-3001
325 cc	12,3 cm	3,8 cm	354-3251	RSZ-3251
350 cc	12,6 cm	3,8 cm	354-3501	RSZ-3501
375 cc	12,9 cm	3,9 cm	354-3751	RSZ-3751
400 cc	13,2 cm	4,0 cm	354-4001	RSZ-4001
450 cc	13,7 cm	4,1 cm	354-4501	RSZ-4501
500 cc	14,1 cm	4,2 cm	354-5001	RSZ-5001
550 cc	14,4 cm	4,4 cm	354-5501	RSZ-5501
600 cc	14,7 cm	4,5 cm	354-6001	RSZ-6001
700 cc	15,7 cm	4,8 cm	354-7001	RSZ-7001
800 cc	16,6 cm	5,0 cm	354-8001	RSZ-8001

SILTEX™ Round Breast Implants,  
Ultra High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
135 cc	8,0 cm	4,3 cm	354-5135	RSZ-5135
160 cc	8,4 cm	4,5 cm	354-5160	RSZ-5160
185 cc	8,6 cm	4,6 cm	354-5185	RSZ-5185
215 cc	8,9 cm	4,7 cm	354-5215	RSZ-5215
240 cc	9,2 cm	4,9 cm	354-5240	RSZ-5240
270 cc	9,5 cm	5,0 cm	354-5270	RSZ-5270
295 cc	9,9 cm	5,2 cm	354-5295	RSZ-5295
320 cc	10,0 cm	5,3 cm	354-5320	RSZ-5320
350 cc	10,3 cm	5,4 cm	354-5350	RSZ-5350
375 cc	10,5 cm	5,5 cm	354-5375	RSZ-5375
400 cc	10,8 cm	5,6 cm	354-5400	RSZ-5400
430 cc	11,1 cm	5,8 cm	354-5430	RSZ-5430
455 cc	11,3 cm	5,9 cm	354-5455	RSZ-5455
480 cc	11,6 cm	6,0 cm	354-5480	RSZ-5480
535 cc	12,1 cm	6,3 cm	354-5535	RSZ-5535
590 cc	12,6 cm	6,5 cm	354-5590	RSZ-5590

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

\*Resterilizable up to 10 times.

Round Gel Breast Implants,  
Cohesive I™ CONT'D

Smooth Round Breast Implants,  
Low Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
100 cc	9,4 cm	1,9 cm	310JSP1100	RSZ-1100
125 cc	10,1 cm	2,0 cm	310JSP1125	RSZ-1125
150 cc	10,7 cm	2,2 cm	310JSP1150	RSZ-1150
175 cc	11,3 cm	2,3 cm	310JSP1175	RSZ-1175
200 cc	11,8 cm	2,4 cm	310JSP1200	RSZ-1200
225 cc	12,3 cm	2,5 cm	310JSP1225	RSZ-1225
250 cc	12,7 cm	2,6 cm	310JSP1250	RSZ-1250
275 cc	13,1 cm	2,7 cm	310JSP1275	RSZ-1275
300 cc	13,5 cm	2,7 cm	310JSP1300	RSZ-1300
325 cc	13,9 cm	2,8 cm	310JSP1325	RSZ-1325
350 cc	14,2 cm	2,9 cm	310JSP1350	RSZ-1350
375 cc	14,5 cm	2,9 cm	310JSP1375	RSZ-1375
400 cc	14,8 cm	3,0 cm	310JSP1400	RSZ-1400
450 cc	15,4 cm	3,1 cm	310JSP1450	RSZ-1450
500 cc	16,0 cm	3,3 cm	310JSP1500	RSZ-1500
550 cc	16,5 cm	3,4 cm	310JSP1550	RSZ-1550
600 cc	17,0 cm	3,5 cm	310JSP1600	RSZ-1600

Smooth Round Breast Implants,  
Moderate Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
100 cc	9,3 cm	2,1 cm	350-7100 BC	RSZ-7100
125 cc	10,0 cm	2,2 cm	350-7125 BC	RSZ-7125
150 cc	10,3 cm	2,3 cm	350-7150 BC	RSZ-7150
175 cc	11,2 cm	2,4 cm	350-7175 BC	RSZ-7175
200 cc	11,7 cm	2,5 cm	350-7200 BC	RSZ-7200
225 cc	12,2 cm	2,6 cm	350-7225 BC	RSZ-7225
250 cc	12,3 cm	2,8 cm	350-7250 BC	RSZ-7250
275 cc	13,2 cm	2,9 cm	350-7275 BC	RSZ-7275
300 cc	13,5 cm	3,0 cm	350-7300 BC	RSZ-7300
325 cc	13,9 cm	3,0 cm	350-7325 BC	RSZ-7325
350 cc	14,2 cm	3,1 cm	350-7350 BC	RSZ-7350
375 cc	14,4 cm	3,2 cm	350-7375 BC	RSZ-7375
400 cc	14,5 cm	3,2 cm	350-7400 BC	RSZ-7400
450 cc	14,9 cm	3,4 cm	350-7450 BC	RSZ-7450
500 cc	15,2 cm	3,6 cm	350-7500 BC	RSZ-7500
550 cc	15,9 cm	3,6 cm	350-7550 BC	RSZ-7550
600 cc	16,5 cm	3,7 cm	350-7600 BC	RSZ-7600
700 cc	17,4 cm	3,9 cm	350-7700 BC	RSZ-7700
800 cc	18,2 cm	4,1 cm	350-7800 BC	RSZ-7800


\*Restorable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Round Gel Breast Implants,  
Cohesive I™ CONT'D

Smooth Round Breast Implants,  
Moderate Plus Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
100 cc	8,2 cm	2,7 cm	350-1001 BC	RSZ-1001
125 cc	8,9 cm	2,8 cm	350-1251 BC	RSZ-1251
150 cc	9,5 cm	2,9 cm	350-1501 BC	RSZ-1501
175 cc	10,0 cm	3,1 cm	350-1751 BC	RSZ-1751
200 cc	10,5 cm	3,2 cm	350-2001 BC	RSZ-2001
225 cc	10,9 cm	3,3 cm	350-2251 BC	RSZ-2251
250 cc	11,3 cm	3,4 cm	350-2501 BC	RSZ-2501
275 cc	11,7 cm	3,5 cm	350-2751 BC	RSZ-2751
300 cc	12,0 cm	3,6 cm	350-3001 BC	RSZ-3001
325 cc	12,3 cm	3,8 cm	350-3251 BC	RSZ-3251
350 cc	12,5 cm	3,9 cm	350-3501 BC	RSZ-3501
375 cc	12,8 cm	4,0 cm	350-3751 BC	RSZ-3751
400 cc	13,1 cm	4,0 cm	350-4001 BC	RSZ-4001
450 cc	13,6 cm	4,2 cm	350-4501 BC	RSZ-4501
500 cc	14,1 cm	4,3 cm	350-5001 BC	RSZ-5001
550 cc	14,6 cm	4,5 cm	350-5501 BC	RSZ-5501
600 cc	15,0 cm	4,6 cm	350-6001 BC	RSZ-6001
700 cc	15,8 cm	4,9 cm	350-7001 BC	RSZ-7001
800 cc	16,5 cm	5,1 cm	350-8001 BC	RSZ-8001

Smooth Round Breast Implants,  
High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
125 cc	8,3 cm	3,5 cm	350-1254 BC	RSZ-1254
150 cc	8,8 cm	3,7 cm	350-1504 BC	RSZ-1504
175 cc	9,3 cm	3,9 cm	350-1754 BC	RSZ-1754
200 cc	9,7 cm	4,0 cm	350-2004 BC	RSZ-2004
225 cc	10,1 cm	4,2 cm	350-2254 BC	RSZ-2254
250 cc	10,5 cm	4,3 cm	350-2504 BC	RSZ-2504
275 cc	10,8 cm	4,4 cm	350-2754 BC	RSZ-2754
300 cc	11,1 cm	4,5 cm	350-3004 BC	RSZ-3004
325 cc	11,4 cm	4,6 cm	350-3254 BC	RSZ-3254
350 cc	11,7 cm	4,8 cm	350-3504 BC	RSZ-3504
375 cc	12,0 cm	4,8 cm	350-3754 BC	RSZ-3754
400 cc	12,2 cm	5,0 cm	350-4004 BC	RSZ-4004
425 cc	12,5 cm	5,0 cm	350-4254 BC	RSZ-4254
450 cc	12,8 cm	5,1 cm	350-4504 BC	RSZ-4504
500 cc	13,2 cm	5,3 cm	350-5004 BC	RSZ-5004
550 cc	13,6 cm	5,5 cm	350-5504 BC	RSZ-5504
600 cc	14,0 cm	5,6 cm	350-6004 BC	RSZ-6004
650 cc	14,4 cm	5,7 cm	350-6504 BC	RSZ-6504
700 cc	14,8 cm	5,8 cm	350-7004 BC	RSZ-7004
800 cc	15,5 cm	6,0 cm	350-8004 BC	RSZ-8004

\*Restorable up to 10 times.

Round Gel Breast Implants,  
Cohesive I™ CONT'D

Smooth Round Breast Implants,  
Ultra High Profile, Cohesive I™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
135 cc	7,8 cm	4,1 cm	350-5135 BC	RSZ-5135
160 cc	8,2 cm	4,3 cm	350-5160 BC	RSZ-5160
185 cc	8,4 cm	4,4 cm	350-5185 BC	RSZ-5185
215 cc	8,7 cm	4,5 cm	350-5215 BC	RSZ-5215
240 cc	9,0 cm	4,7 cm	350-5240 BC	RSZ-5240
270 cc	9,3 cm	4,8 cm	350-5270 BC	RSZ-5270
295 cc	9,7 cm	5,0 cm	350-5295 BC	RSZ-5295
320 cc	9,8 cm	5,1 cm	350-5320 BC	RSZ-5320
350 cc	10,1 cm	5,2 cm	350-5350 BC	RSZ-5350
375 cc	10,4 cm	5,3 cm	350-5375 BC	RSZ-5375
400 cc	10,6 cm	5,4 cm	350-5400 BC	RSZ-5400
430 cc	10,9 cm	5,6 cm	350-5430 BC	RSZ-5430
455 cc	11,2 cm	5,7 cm	350-5455 BC	RSZ-5455
480 cc	11,4 cm	5,8 cm	350-5480 BC	RSZ-5480
535 cc	12,0 cm	6,1 cm	350-5535 BC	RSZ-5535
590 cc	12,5 cm	6,3 cm	350-5590 BC	RSZ-5590

Round Gel Breast Implants,  
Cohesive II™

SILTEX™ Round Breast Implants,  
Moderate Profile, Cohesive II™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
100 cc	9,1 cm	2,6 cm	354-1009	RSZ-7100
125 cc	9,3 cm	2,8 cm	354-1259	RSZ-7125
150 cc	10,1 cm	2,4 cm	354-1509	RSZ-7150
175 cc	10,5 cm	2,9 cm	354-1759	RSZ-7175
200 cc	11,2 cm	2,8 cm	354-2009	RSZ-7200
225 cc	11,5 cm	3,0 cm	354-2259	RSZ-7225
250 cc	11,6 cm	3,5 cm	354-2509	RSZ-7250
275 cc	12,6 cm	3,1 cm	354-2759	RSZ-7275
300 cc	12,6 cm	3,6 cm	354-3009	RSZ-7300
325 cc	13,1 cm	3,6 cm	354-3259	RSZ-7325
350 cc	13,4 cm	3,5 cm	354-3509	RSZ-7350
375 cc	13,4 cm	3,9 cm	354-3759	RSZ-7375
400 cc	13,7 cm	3,8 cm	354-4009	RSZ-7400
450 cc	13,7 cm	4,6 cm	354-4509	RSZ-7450
500 cc	14,2 cm	4,8 cm	354-5009	RSZ-7500
550 cc	14,5 cm	4,5 cm	354-5509	RSZ-7550
600 cc	15,3 cm	4,6 cm	354-6009	RSZ-7600
700 cc	16,3 cm	5,0 cm	354-7009	RSZ-7700
800 cc	16,7 cm	5,4 cm	354-8009	RSZ-7800

\*Restorable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Round Gel Breast Implants,  
Cohesive II™ CONT'D

SILTEX™ Round Breast Implants,  
Moderate Plus Profile, Cohesive II™



Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
100 cc	8,1 cm	2,7 cm	324-5100	RSZ-1001
125 cc	8,8 cm	2,9 cm	324-5125	RSZ-1251
150 cc	9,4 cm	3,0 cm	324-5150	RSZ-1501
175 cc	10,0 cm	3,2 cm	324-5175	RSZ-1751
200 cc	10,5 cm	3,3 cm	324-5200	RSZ-2001
225 cc	10,9 cm	3,5 cm	324-5225	RSZ-2251
250 cc	11,3 cm	3,6 cm	324-5250	RSZ-2501
275 cc	11,7 cm	3,7 cm	324-5275	RSZ-2751
300 cc	12,0 cm	3,7 cm	324-5300	RSZ-3001
325 cc	12,3 cm	3,8 cm	324-5325	RSZ-3251
350 cc	12,6 cm	3,8 cm	324-5350	RSZ-3501
375 cc	12,9 cm	3,9 cm	324-5375	RSZ-3751
400 cc	13,2 cm	4,0 cm	324-5400	RSZ-4001
450 cc	13,7 cm	4,1 cm	324-5450	RSZ-4501
500 cc	14,1 cm	4,2 cm	324-5500	RSZ-5001
550 cc	14,4 cm	4,4 cm	324-5550	RSZ-5501
600 cc	14,7 cm	4,5 cm	324-5600	RSZ-6001
700 cc	15,7 cm	4,8 cm	324-5700	RSZ-7001
800 cc	16,6 cm	5,0 cm	324-5800	RSZ-8001

SILTEX™ Round Breast Implants,  
High Profile, Cohesive II™




Volume	Diameter	Projection	Catalog #	Restorable Gel Sizer*
125 cc	8,4 cm	3,6 cm	324-4125	RSZ-1254
150 cc	8,9 cm	3,8 cm	324-4150	RSZ-1504
175 cc	9,4 cm	4,0 cm	324-4175	RSZ-1754
200 cc	9,9 cm	4,1 cm	324-4200	RSZ-2004
225 cc	10,2 cm	4,3 cm	324-4225	RSZ-2254
250 cc	10,5 cm	4,5 cm	324-4250	RSZ-2504
275 cc	10,9 cm	4,6 cm	324-4275	RSZ-2754
300 cc	11,1 cm	4,7 cm	324-4300	RSZ-3004
325 cc	11,5 cm	4,8 cm	324-4325	RSZ-3254
350 cc	11,7 cm	4,9 cm	324-4350	RSZ-3504
375 cc	12,0 cm	5,0 cm	324-4375	RSZ-3754
400 cc	12,3 cm	5,1 cm	324-4400	RSZ-4004
425 cc	12,5 cm	5,2 cm	324-4425	RSZ-4254
450 cc	12,7 cm	5,2 cm	324-4450	RSZ-4504
500 cc	13,2 cm	5,4 cm	324-4500	RSZ-5004
550 cc	13,5 cm	5,6 cm	324-4550	RSZ-5504
600 cc	14,0 cm	5,7 cm	324-4600	RSZ-6004
650 cc	14,3 cm	5,8 cm	324-4650	RSZ-6504
700 cc	14,7 cm	6,0 cm	324-4700	RSZ-7004
800 cc	15,4 cm	6,3 cm	324-4800	RSZ-8004

\*Restorable up to 10 times.



CPG™ BREAST IMPLANT NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate Plus 2
	Low 1	Moderate 1

KEY: DIGIT 1 = COHESIVE LEVEL  
 DIGIT 2 = HEIGHT OF IMPLANT  
 DIGIT 3 = PROJECTION OF IMPLANT



### CPG™ Implants 331, Cohesive III™, Tall Height, Moderate Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	9,2 cm	3,2 cm	2,4 cm	3,9 cm	334-0903	RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	2,6 cm	4,1 cm	334-0953	RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	2,7 cm	4,3 cm	334-1003	RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	2,8 cm	4,6 cm	334-1053	RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	3,0 cm	4,8 cm	334-1103	RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	3,1 cm	5,0 cm	334-1153	RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	3,2 cm	5,2 cm	334-1203	RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	3,4 cm	5,4 cm	334-1253	RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	3,5 cm	5,6 cm	334-1303	RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	3,6 cm	5,9 cm	334-1353	RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	3,8 cm	6,1 cm	334-1403	RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	3,9 cm	6,3 cm	334-1453	RSZ-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	4,0 cm	6,5 cm	334-1503	RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	4,2 cm	6,7 cm	334-1553	RSZ-1553

### CPG™ Implants 321, Cohesive III™, Medium Height, Moderate Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
120 cc	9,0 cm	8,5 cm	3,3 cm	3,1 cm	4,3 cm	354-0908	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	3,2 cm	4,5 cm	354-0958	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	3,4 cm	4,7 cm	354-1008	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	3,6 cm	5,0 cm	354-1058	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	3,6 cm	5,1 cm	354-1108	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	3,8 cm	5,3 cm	354-1158	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	4,0 cm	5,6 cm	354-1208	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	4,1 cm	5,8 cm	354-1258	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	4,3 cm	6,0 cm	354-1308	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	4,5 cm	6,3 cm	354-1358	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	4,6 cm	6,5 cm	354-1408	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	5,0 cm	6,9 cm	354-1458	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	5,1 cm	7,2 cm	354-1508	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	5,5 cm	7,6 cm	354-1608	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	5,8 cm	8,1 cm	354-1708	RSZ-1708

### CPG™ Implants 311, Cohesive III™, Low Height, Moderate Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
100 cc	9,0 cm	8,0 cm	3,0 cm	2,5 cm	4,1 cm	334-0901	RSZ-0090
120 cc	9,5 cm	8,5 cm	3,2 cm	2,6 cm	4,3 cm	334-0951	RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	2,8 cm	4,5 cm	334-1001	RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	2,9 cm	4,7 cm	334-1051	RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	3,0 cm	5,0 cm	334-1101	RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	3,2 cm	5,2 cm	334-1151	RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	3,3 cm	5,4 cm	334-1201	RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	3,5 cm	5,6 cm	334-1251	RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	3,6 cm	5,9 cm	334-1301	RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	3,7 cm	6,1 cm	334-1351	RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	3,9 cm	6,3 cm	334-1401	RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	4,0 cm	6,5 cm	334-1451	RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	4,1 cm	6,8 cm	334-1501	RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	4,3 cm	7,0 cm	334-1551	RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	4,4 cm	7,2 cm	334-1601	RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	4,6 cm	7,4 cm	334-1651	RSZ-0165

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

CPG™ Gel Breast Implants, Cohesive III™ CONT'D

CPG™ Gel Breast Implants, Cohesive III™ CONT'D

CPG™ Implants 332,  
Cohesive III™, Tall Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
145 cc	9,0 cm	9,4 cm	3,8 cm	2,8 cm	4,6 cm	334-0909	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	3,0 cm	4,9 cm	334-0959	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	3,1 cm	5,1 cm	334-1009	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	3,3 cm	5,4 cm	334-1059	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	3,5 cm	5,7 cm	334-1109	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	3,6 cm	5,9 cm	334-1159	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	3,8 cm	6,2 cm	334-1209	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	3,9 cm	6,4 cm	334-1259	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	4,1 cm	6,7 cm	334-1309	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	4,2 cm	7,0 cm	334-1359	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	4,4 cm	7,2 cm	334-1409	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	4,6 cm	7,5 cm	334-1459	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	4,7 cm	7,7 cm	334-1509	RSZ-1509

CPG™ Implants 333,  
Cohesive III™, Tall Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
180 cc	9,0 cm	9,4 cm	4,6 cm	3,1 cm	5,5 cm	334-0904	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	3,3 cm	5,8 cm	334-0954	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	3,4 cm	6,1 cm	334-1004	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	3,6 cm	6,4 cm	334-1054	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	3,8 cm	6,8 cm	334-1104	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	3,9 cm	7,1 cm	334-1154	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	4,1 cm	7,4 cm	334-1204	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	4,3 cm	7,7 cm	334-1254	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	4,5 cm	8,0 cm	334-1304	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	4,6 cm	8,3 cm	334-1354	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	4,8 cm	8,6 cm	334-1404	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	5,0 cm	8,9 cm	334-1454	RSZ-4145

CPG™ Implants 322,  
Cohesive III™, Medium Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
140 cc	9,0 cm	8,5 cm	3,8 cm	2,9 cm	4,7 cm	334-0905	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	3,0 cm	4,9 cm	334-0955	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	3,2 cm	5,2 cm	334-1005	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	3,3 cm	5,5 cm	334-1055	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	3,5 cm	5,7 cm	334-1105	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	3,7 cm	6,0 cm	334-1155	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	3,8 cm	6,3 cm	334-1205	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	4,0 cm	6,5 cm	334-1255	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	4,1 cm	6,8 cm	334-1305	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	4,3 cm	7,0 cm	334-1355	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	4,5 cm	7,3 cm	334-1405	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	4,6 cm	7,6 cm	334-1455	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	4,8 cm	7,8 cm	334-1505	RSZ-1505

CPG™ Implants 323,  
Cohesive III™, Medium Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
165 cc	9,0 cm	8,5 cm	4,6 cm	3,0 cm	5,3 cm	334-0902	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	3,2 cm	5,6 cm	334-0952	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	3,4 cm	5,9 cm	334-1002	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	3,5 cm	6,2 cm	334-1052	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	3,7 cm	6,5 cm	334-1102	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	3,9 cm	6,8 cm	334-1152	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	4,0 cm	7,1 cm	334-1202	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	4,2 cm	7,4 cm	334-1252	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	4,4 cm	7,7 cm	334-1302	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	4,5 cm	8,0 cm	334-1352	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	4,7 cm	8,3 cm	334-1402	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	4,9 cm	8,6 cm	334-1452	RSZ-1452

CPG™ Implants 312,  
Cohesive III™, Low Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	8,0 cm	3,8 cm	2,5 cm	4,6 cm	334-0907	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	2,6 cm	4,8 cm	334-0957	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	2,7 cm	5,1 cm	334-1007	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	2,9 cm	5,4 cm	334-1057	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	3,0 cm	5,6 cm	334-1107	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	3,1 cm	5,9 cm	334-1157	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	3,3 cm	6,1 cm	334-1207	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	3,4 cm	6,4 cm	334-1257	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	3,6 cm	6,6 cm	334-1307	RSZ-1307
415 cc	13,5 cm	12,0 cm	5,7 cm	3,7 cm	6,9 cm	334-1357	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	3,8 cm	7,1 cm	334-1407	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	4,0 cm	7,4 cm	334-1457	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	4,1 cm	7,7 cm	334-1507	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	4,4 cm	8,2 cm	334-1607	RSZ-1607

CPG™ Implants 313,  
Cohesive III™, Low Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #	Resterilizable Gel Sizer*
130 cc	9,0 cm	8,1 cm	4,4 cm	2,6 cm	5,2 cm	334-0906	RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	2,7 cm	5,4 cm	334-0956	RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	2,9 cm	5,7 cm	334-1006	RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	3,0 cm	6,0 cm	334-1056	RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	3,2 cm	6,3 cm	334-1106	RSZ-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	3,3 cm	6,6 cm	334-1156	RSZ-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	3,5 cm	6,9 cm	334-1206	RSZ-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	3,6 cm	7,2 cm	334-1256	RSZ-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	3,8 cm	7,4 cm	334-1306	RSZ-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	3,9 cm	7,7 cm	334-1356	RSZ-1356
490 cc	14,0 cm	12,6 cm	6,1 cm	4,0 cm	8,0 cm	334-1406	RSZ-1406
545 cc	14,5 cm	13,1 cm	6,3 cm	4,2 cm	8,3 cm	334-1456	RSZ-1456
605 cc	15,0 cm	13,5 cm	6,6 cm	4,3 cm	8,6 cm	334-1506	RSZ-1506

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

\*Resterilizable up to 10 times.



BECKER™ Expander/Breast Implants

SILTEX™ Round BECKER™ 25  
Expander/Implants, Cohesive I™


(25% Silicone Gel, Cohesive I™, in outer lumen, 75% Saline in inner lumen)



Volume	Gel Volume	TEMPORARY OVEREXPANSION VOLUMES		Total Saline	Total Gel-Saline	Diameter	Proj.	Cat.#
		Max Saline	Total Gel-Saline					
150 cc	40 cc	185 cc	225 cc	85-150 cc	125-190 cc	10,7 cm	2,4 cm	354-1500
200 cc	50 cc	250 cc	300 cc	125-200 cc	175-250 cc	11,6 cm	2,7 cm	354-2000
250 cc	60 cc	315 cc	375 cc	165-255 cc	225-315 cc	12,1 cm	3,0 cm	354-2500
300 cc	75 cc	375 cc	450 cc	200-300 cc	275-375 cc	12,9 cm	3,0 cm	354-3000
350 cc	90 cc	435 cc	525 cc	235-350 cc	325-440 cc	13,5 cm	3,0 cm	354-3500
400 cc	100 cc	500 cc	600 cc	275-400 cc	375-500 cc	13,7 cm	3,4 cm	354-4000
500 cc	125 cc	625 cc	750 cc	350-500 cc	475-625 cc	14,7 cm	3,7 cm	354-5000
600 cc	150 cc	750 cc	900 cc	425-600 cc	575-750 cc	15,6 cm	4,2 cm	354-6000
700 cc	175 cc	875 cc	1050 cc	500-700 cc	675-875 cc	16,8 cm	4,3 cm	354-7000
800 cc	200 cc	1000 cc	1200 cc	575-800 cc	775-1000 cc	17,2 cm	4,6 cm	354-8000

SILTEX™ Round BECKER™ 50  
Expander/Implants, Cohesive I™

(50% Silicone Gel, Cohesive I™, in outer lumen, 50% Saline in inner lumen)



Volume	Gel Volume	Total Saline	Total Gel-Saline	Diameter	Proj.	Catalog #
300 cc	150 cc	150-200 cc	300-350 cc	12,9 cm	3,0 cm	354-1515
400 cc	200 cc	200-300 cc	400-500 cc	13,9 cm	3,4 cm	354-2020
500 cc	250 cc	250-350 cc	500-600 cc	14,7 cm	3,7 cm	354-2525
600 cc	300 cc	300-425 cc	600-725 cc	15,6 cm	4,2 cm	354-3030
700 cc	350 cc	350-500 cc	700-850 cc	16,0 cm	5,0 cm	354-3535

SILTEX™ CONTOUR PROFILE™  
BECKER™ 35 Expander/Implants, Cohesive II™

(35% Silicone Gel, Cohesive II™, in outer lumen, 65% Saline in inner lumen)



Volume	Gel Volume	TEMPORARY OVEREXPANSION VOLUMES		Total Saline	Total Gel-Saline	Width	Height	Proj.	Cat.#
		Max Saline	Total Gel-Saline						
145 cc	50 cc	120 cc	170 cc	85-95 cc	135-145 cc	9,5 cm	8,3 cm	4,2 cm	324-0955
195 cc	70 cc	155 cc	225 cc	110-125 cc	180-195 cc	10,5 cm	9,1 cm	4,6 cm	324-1055
255 cc	90 cc	205 cc	295 cc	145-165 cc	235-255 cc	11,5 cm	10,0 cm	5,1 cm	324-1155
290 cc	100 cc	235 cc	335 cc	170-190 cc	270-290 cc	12,0 cm	10,5 cm	5,3 cm	324-1205
325 cc	110 cc	270 cc	380 cc	190-215 cc	300-325 cc	12,5 cm	10,9 cm	5,5 cm	324-1255
365 cc	125 cc	300 cc	425 cc	215-240 cc	340-365 cc	13,0 cm	11,3 cm	5,7 cm	324-1305
400 cc	130 cc	335 cc	465 cc	240-270 cc	370-400 cc	13,5 cm	11,8 cm	5,9 cm	324-1355
460 cc	160 cc	375 cc	535 cc	270-300 cc	430-460 cc	14,0 cm	12,2 cm	6,2 cm	324-1405
565 cc	195 cc	460 cc	655 cc	330-370 cc	525-565 cc	15,0 cm	13,1 cm	6,6 cm	324-1505
685 cc	240 cc	555 cc	795 cc	400-445 cc	640-685 cc	16,0 cm	13,9 cm	7,0 cm	324-1605

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

SPECTRA™ Adjustable Gel Breast Implants

SILTEX™ Round SPECTRA™  
Adjustable Gel Breast Implants, Cohesive I™‡



Volume	Gel Volume	Saline Volume	NOMINAL FILL		MAXIMUM FILL		Catalog #
			Diameter	Projection	Diameter	Projection	
190 cc	180 cc	10-50cc	9,9 cm	4,1 cm	9,4 cm	5,1 cm	314-4190
215 cc	205 cc	10-50cc	10,2 cm	4,3 cm	9,8 cm	5,2 cm	314-4215
240 cc	230 cc	10-50cc	10,5 cm	4,5 cm	10,2 cm	5,3 cm	314-4240
260 cc	250 cc	10-50cc	10,9 cm	4,6 cm	10,5 cm	5,4 cm	314-4260
285 cc	275 cc	10-50cc	11,1 cm	4,7 cm	10,8 cm	5,5 cm	314-4285
310 cc	300 cc	10-50cc	11,5 cm	4,8 cm	11,1 cm	5,6 cm	314-4310
335 cc	325 cc	10-50cc	11,7 cm	4,9 cm	11,4 cm	5,7 cm	314-4335
355 cc	345 cc	10-50cc	12,0 cm	5,0 cm	11,7 cm	5,8 cm	314-4355
380 cc	370 cc	10-50cc	12,3 cm	5,1 cm	11,9 cm	5,8 cm	314-4380
405 cc	395 cc	10-50cc	12,5 cm	5,2 cm	12,2 cm	5,9 cm	314-4405
430 cc	420 cc	10-50cc	12,7 cm	5,2 cm	12,4 cm	6,0 cm	314-4430
480 cc	465 cc	15-75cc	13,2 cm	5,4 cm	12,9 cm	6,1 cm	314-4480
525 cc	510 cc	15-75cc	13,5 cm	5,6 cm	13,3 cm	6,2 cm	314-4525
575 cc	560 cc	15-75cc	14,0 cm	5,7 cm	13,8 cm	6,4 cm	314-4575
620 cc	605 cc	15-100cc	14,3 cm	5,8 cm	14,1 cm	6,5 cm	314-4620
665 cc	650 cc	15-100cc	14,7 cm	6,0 cm	14,5 cm	6,6 cm	314-4665
760 cc	745 cc	15-100cc	15,4 cm	6,3 cm	15,1 cm	6,7 cm	314-4760

Smooth Round SPECTRA™  
Adjustable Gel Breast Implants, Cohesive I™‡



Volume	Gel Volume	Saline Volume	NOMINAL FILL		MAXIMUM FILL		Catalog #
			Diameter	Projection	Diameter	Projection	
190 cc	180 cc	10-50cc	9,7 cm	4,0 cm	9,4 cm	5,1 cm	310-4190
215 cc	205 cc	10-50cc	10,1 cm	4,2 cm	9,8 cm	5,2 cm	310-4215
240 cc	230 cc	10-50cc	10,5 cm	4,3 cm	10,2 cm	5,3 cm	310-4240
260 cc	250 cc	10-50cc	10,8 cm	4,4 cm	10,5 cm	5,4 cm	310-4260
285 cc	275 cc	10-50cc	11,1 cm	4,5 cm	10,8 cm	5,5 cm	310-4285
310 cc	300 cc	10-50cc	11,4 cm	4,6 cm	11,1 cm	5,6 cm	310-4310
335 cc	325 cc	10-50cc	11,7 cm	4,8 cm	11,4 cm	5,7 cm	310-4335
355 cc	345 cc	10-50cc	12,0 cm	4,8 cm	11,7 cm	5,8 cm	310-4355
380 cc	370 cc	10-50cc	12,2 cm	5,0 cm	11,9 cm	5,8 cm	310-4380
405 cc	395 cc	10-50cc	12,5 cm	5,0 cm	12,2 cm	5,9 cm	310-4405
430 cc	420 cc	10-50cc	12,8 cm	5,1 cm	12,4 cm	6,0 cm	310-4430
480 cc	465 cc	15-75cc	13,2 cm	5,3 cm	12,9 cm	6,1 cm	310-4480
525 cc	510 cc	15-75cc	13,6 cm	5,5 cm	13,3 cm	6,2 cm	310-4525
575 cc	560 cc	15-75cc	14,0 cm	5,6 cm	13,8 cm	6,4 cm	310-4575
620 cc	605 cc	15-100cc	14,4 cm	5,7 cm	14,1 cm	6,5 cm	310-4620
665 cc	650 cc	15-100cc	14,8 cm	5,8 cm	14,5 cm	6,6 cm	310-4665
760 cc	745 cc	15-100cc	15,5 cm	6,0 cm	15,1 cm	6,7 cm	310-4760

‡Dome Pack for postoperative adjustment can be purchased separately

Sterile Resterilizable  
Gel Breast Implant Sizers\*

Sterile Resterilizable  
Gel Breast Implant Sizers\* CONT'D

Round Gel Breast Implant Sizers,  
Low Profile



Volume	Diameter	Projection	Catalog #
100 cc	9,4 cm	1,9 cm	RSZ-1100
125 cc	10,1 cm	2,0 cm	RSZ-1125
150 cc	10,7 cm	2,2 cm	RSZ-1150
175 cc	11,3 cm	2,3 cm	RSZ-1175
200 cc	11,8 cm	2,4 cm	RSZ-1200
225 cc	12,3 cm	2,5 cm	RSZ-1225
250 cc	12,7 cm	2,6 cm	RSZ-1250
275 cc	13,1 cm	2,7 cm	RSZ-1275
300 cc	13,5 cm	2,7 cm	RSZ-1300
325 cc	13,9 cm	2,8 cm	RSZ-1325
350 cc	14,2 cm	2,9 cm	RSZ-1350
375 cc	14,5 cm	2,9 cm	RSZ-1375
400 cc	14,8 cm	3,0 cm	RSZ-1400
450 cc	15,4 cm	3,1 cm	RSZ-1450
500 cc	16,0 cm	3,3 cm	RSZ-1500
550 cc	16,5 cm	3,4 cm	RSZ-1550
600 cc	17,0 cm	3,5 cm	RSZ-1600

Round Gel Breast Implant Sizers,  
Moderate Plus Profile



Volume	Diameter	Projection	Catalog #
100 cc	8,2 cm	2,7 cm	RSZ-1001
125 cc	8,9 cm	2,8 cm	RSZ-1251
150 cc	9,5 cm	2,9 cm	RSZ-1501
175 cc	10,0 cm	3,1 cm	RSZ-1751
200 cc	10,5 cm	3,2 cm	RSZ-2001
225 cc	10,9 cm	3,3 cm	RSZ-2251
250 cc	11,3 cm	3,4 cm	RSZ-2501
275 cc	11,7 cm	3,5 cm	RSZ-2751
300 cc	12,0 cm	3,6 cm	RSZ-3001
325 cc	12,3 cm	3,8 cm	RSZ-3251
350 cc	12,5 cm	3,9 cm	RSZ-3501
375 cc	12,8 cm	4,0 cm	RSZ-3751
400 cc	13,1 cm	4,0 cm	RSZ-4001
450 cc	13,6 cm	4,2 cm	RSZ-4501
500 cc	14,1 cm	4,3 cm	RSZ-5001
550 cc	14,6 cm	4,5 cm	RSZ-5501
600 cc	15,0 cm	4,6 cm	RSZ-6001
700 cc	15,8 cm	4,9 cm	RSZ-7001
800 cc	16,5 cm	5,1 cm	RSZ-8001

Round Gel Breast Implant Sizers,  
Moderate Profile



Volume	Diameter	Projection	Catalog #
100 cc	9,3 cm	2,1 cm	RSZ-7100
125 cc	10,0 cm	2,2 cm	RSZ-7125
150 cc	10,3 cm	2,3 cm	RSZ-7150
175 cc	11,2 cm	2,4 cm	RSZ-7175
200 cc	11,7 cm	2,5 cm	RSZ-7200
225 cc	12,2 cm	2,6 cm	RSZ-7225
250 cc	12,3 cm	2,8 cm	RSZ-7250
275 cc	13,2 cm	2,9 cm	RSZ-7275
300 cc	13,5 cm	3,0 cm	RSZ-7300
325 cc	13,9 cm	3,0 cm	RSZ-7325
350 cc	14,2 cm	3,1 cm	RSZ-7350
375 cc	14,4 cm	3,2 cm	RSZ-7375
400 cc	14,5 cm	3,2 cm	RSZ-7400
450 cc	14,9 cm	3,4 cm	RSZ-7450
500 cc	15,2 cm	3,6 cm	RSZ-7500
550 cc	15,9 cm	3,6 cm	RSZ-7550
600 cc	16,5 cm	3,7 cm	RSZ-7600
700 cc	17,4 cm	3,9 cm	RSZ-7700
800 cc	18,2 cm	4,1 cm	RSZ-7800

Round Gel Breast Implant Sizers,  
High Profile



Volume	Diameter	Projection	Catalog #
125 cc	8,3 cm	3,5 cm	RSZ-1254
150 cc	8,8 cm	3,7 cm	RSZ-1504
175 cc	9,3 cm	3,9 cm	RSZ-1754
200 cc	9,7 cm	4,0 cm	RSZ-2004
225 cc	10,1 cm	4,2 cm	RSZ-2254
250 cc	10,5 cm	4,3 cm	RSZ-2504
275 cc	10,8 cm	4,4 cm	RSZ-2754
300 cc	11,1 cm	4,5 cm	RSZ-3004
325 cc	11,4 cm	4,6 cm	RSZ-3254
350 cc	11,7 cm	4,8 cm	RSZ-3504
375 cc	12,0 cm	4,8 cm	RSZ-3754
400 cc	12,2 cm	5,0 cm	RSZ-4004
425 cc	12,5 cm	5,0 cm	RSZ-4254
450 cc	12,8 cm	5,1 cm	RSZ-4504
500 cc	13,2 cm	5,3 cm	RSZ-5004
550 cc	13,6 cm	5,5 cm	RSZ-5504
600 cc	14,0 cm	5,6 cm	RSZ-6004
650 cc	14,4 cm	5,7 cm	RSZ-6504
700 cc	14,8 cm	5,8 cm	RSZ-7004
800 cc	15,5 cm	6,0 cm	RSZ-8004

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

\*Resterilizable up to 10 times.

Sterile Resterilizable  
Gel Breast Implant Sizers\* CONT'D

Round Gel Breast Implant Sizers,  
Ultra High Profile



Volume	Diameter	Projection	Catalog #
135 cc	7,8 cm	4,1 cm	RSZ-5135
160 cc	8,2 cm	4,3 cm	RSZ-5160
185 cc	8,4 cm	4,4 cm	RSZ-5185
215 cc	8,7 cm	4,5 cm	RSZ-5215
240 cc	9,0 cm	4,7 cm	RSZ-5240
270 cc	9,3 cm	4,8 cm	RSZ-5270
295 cc	9,7 cm	5,0 cm	RSZ-5295
320 cc	9,8 cm	5,1 cm	RSZ-5320
350 cc	10,1 cm	5,2 cm	RSZ-5350
375 cc	10,4 cm	5,3 cm	RSZ-5375
400 cc	10,6 cm	5,4 cm	RSZ-5400
430 cc	10,9 cm	5,6 cm	RSZ-5430
455 cc	11,2 cm	5,7 cm	RSZ-5455
480 cc	11,4 cm	5,8 cm	RSZ-5480
535 cc	12,0 cm	6,1 cm	RSZ-5535
590 cc	12,5 cm	6,3 cm	RSZ-5590

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Sterile Resterilizable  
Gel Breast Implant Sizers\* CONT'D

CPG™ GEL BREAST IMPLANT SIZERS



CPG™ Implant Sizers 331  
Cohesive III™, Tall Height, Moderate Projection

Volume	Width	Height	Projection	A	Arc Length	Catalog #
125 cc	9,0 cm	9,2 cm	3,2 cm	2,4 cm	3,9 cm	RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	2,6 cm	4,1 cm	RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	2,7 cm	4,3 cm	RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	2,8 cm	4,6 cm	RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	3,0 cm	4,8 cm	RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	3,1 cm	5,0 cm	RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	3,2 cm	5,2 cm	RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	3,4 cm	5,4 cm	RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	3,5 cm	5,6 cm	RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	3,6 cm	5,9 cm	RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	3,8 cm	6,1 cm	RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	3,9 cm	6,3 cm	RSZ-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	4,0 cm	6,5 cm	RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	4,2 cm	6,7 cm	RSZ-1553

CPG™ Implant Sizers 321  
Cohesive III™, Medium Height, Moderate Projection

Volume	Width	Height	Projection	A	Arc Length	Catalog #
120 cc	9,0 cm	8,5 cm	3,3 cm	3,1 cm	4,3 cm	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	3,2 cm	4,5 cm	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	3,4 cm	4,7 cm	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	3,6 cm	5,0 cm	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	3,6 cm	5,1 cm	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	3,8 cm	5,3 cm	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	4,0 cm	5,6 cm	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	4,1 cm	5,8 cm	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	4,3 cm	6,0 cm	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	4,5 cm	6,3 cm	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	4,6 cm	6,5 cm	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	5,0 cm	6,9 cm	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	5,1 cm	7,2 cm	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	5,5 cm	7,6 cm	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	5,8 cm	8,1 cm	RSZ-1708

CPG™ Implant Sizers 311  
Cohesive III™, Low Height, Moderate Projection

Volume	Width	Height	Projection	A	Arc Length	Catalog #
100 cc	9,0 cm	8,0 cm	3,0 cm	2,5 cm	4,1 cm	RSZ-0090
120 cc	9,5 cm	8,5 cm	3,2 cm	2,6 cm	4,3 cm	RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	2,8 cm	4,5 cm	RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	2,9 cm	4,7 cm	RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	3,0 cm	5,0 cm	RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	3,2 cm	5,2 cm	RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	3,3 cm	5,4 cm	RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	3,5 cm	5,6 cm	RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	3,6 cm	5,9 cm	RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	3,7 cm	6,1 cm	RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	3,9 cm	6,3 cm	RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	4,0 cm	6,5 cm	RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	4,1 cm	6,8 cm	RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	4,3 cm	7,0 cm	RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	4,4 cm	7,2 cm	RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	4,6 cm	7,4 cm	RSZ-0165

\*Resterilizable up to 10 times.

Sterile Resterilizable  
Gel Breast Implant Sizers\* CONT'D

Sterile Resterilizable  
Gel Breast Implant Sizers\* CONT'D

CPG™ GEL BREAST IMPLANT SIZERS

CPG™ GEL BREAST IMPLANT SIZERS

CPG™ Implant Sizers 332  
Cohesive III™, Tall Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
145 cc	9,0 cm	9,4 cm	3,8 cm	2,8 cm	4,6 cm	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	3,0 cm	4,9 cm	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	3,1 cm	5,1 cm	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	3,3 cm	5,4 cm	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	3,5 cm	5,7 cm	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	3,6 cm	5,9 cm	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	3,8 cm	6,2 cm	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	3,9 cm	6,4 cm	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	4,1 cm	6,7 cm	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	4,2 cm	7,0 cm	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	4,4 cm	7,2 cm	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	4,6 cm	7,5 cm	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	4,7 cm	7,7 cm	RSZ-1509

CPG™ Implant Sizers 333  
Cohesive III™, Tall Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
180 cc	9,0 cm	9,4 cm	4,6 cm	3,1 cm	5,5 cm	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	3,3 cm	5,8 cm	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	3,4 cm	6,1 cm	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	3,6 cm	6,4 cm	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	3,8 cm	6,8 cm	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	3,9 cm	7,1 cm	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	4,1 cm	7,4 cm	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	4,3 cm	7,7 cm	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	4,5 cm	8,0 cm	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	4,6 cm	8,3 cm	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	4,8 cm	8,6 cm	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	5,0 cm	8,9 cm	RSZ-4145

CPG™ Implant Sizers 322  
Cohesive III™, Medium Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
140 cc	9,0 cm	8,5 cm	3,8 cm	2,9 cm	4,7 cm	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	3,0 cm	4,9 cm	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	3,2 cm	5,2 cm	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	3,3 cm	5,5 cm	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	3,5 cm	5,7 cm	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	3,7 cm	6,0 cm	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	3,8 cm	6,3 cm	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	4,0 cm	6,5 cm	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	4,1 cm	6,8 cm	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	4,3 cm	7,0 cm	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	4,5 cm	7,3 cm	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	4,6 cm	7,6 cm	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	4,8 cm	7,8 cm	RSZ-1505

CPG™ Implant Sizers 323  
Cohesive III™, Medium Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
165 cc	9,0 cm	8,5 cm	4,6 cm	3,0 cm	5,3 cm	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	3,2 cm	5,6 cm	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	3,4 cm	5,9 cm	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	3,5 cm	6,2 cm	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	3,7 cm	6,5 cm	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	3,9 cm	6,8 cm	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	4,0 cm	7,1 cm	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	4,2 cm	7,4 cm	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	4,4 cm	7,7 cm	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	4,5 cm	8,0 cm	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	4,7 cm	8,3 cm	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	4,9 cm	8,6 cm	RSZ-1452

CPG™ Implant Sizers 312  
Cohesive III™, Low Height, Moderate Plus Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
125 cc	9,0 cm	8,0 cm	3,8 cm	2,5 cm	4,6 cm	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	2,6 cm	4,8 cm	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	2,7 cm	5,1 cm	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	2,9 cm	5,4 cm	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	3,0 cm	5,6 cm	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	3,1 cm	5,9 cm	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	3,3 cm	6,1 cm	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	3,4 cm	6,4 cm	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	3,6 cm	6,6 cm	RSZ-1307
415 cc	13,5 cm	12,0 cm	5,7 cm	3,7 cm	6,9 cm	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	3,8 cm	7,1 cm	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	4,0 cm	7,4 cm	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	4,1 cm	7,7 cm	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	4,4 cm	8,2 cm	RSZ-1607

CPG™ Implant Sizers 313  
Cohesive III™, Low Height, High Projection



Volume	Width	Height	Projection	A	Arc Length	Catalog #
130 cc	9,0 cm	8,1 cm	4,4 cm	2,6 cm	5,2 cm	RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	2,7 cm	5,4 cm	RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	2,9 cm	5,7 cm	RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	3,0 cm	6,0 cm	RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	3,2 cm	6,3 cm	RSZ-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	3,3 cm	6,6 cm	RSZ-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	3,5 cm	6,9 cm	RSZ-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	3,6 cm	7,2 cm	RSZ-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	3,8 cm	7,4 cm	RSZ-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	3,9 cm	7,7 cm	RSZ-1356
490 cc	14,0 cm	12,6 cm	6,1 cm	4,0 cm	8,0 cm	RSZ-1406
545 cc	14,5 cm	13,1 cm	6,3 cm	4,2 cm	8,3 cm	RSZ-1456
605 cc	15,0 cm	13,5 cm	6,6 cm	4,3 cm	8,6 cm	RSZ-1506

\*Resterilizable up to 10 times.

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

\*Resterilizable up to 10 times.

## CONTOUR PROFILE™ Tissue Expanders

**SILTEX™ CONTOUR PROFILE™**  
 Tissue Expanders, Low Height,  
 Style 6100\*


Volume	Width	Height	Projection	Catalog #
250 cc	11,4 cm	8,1 cm	6,1 cm	354-6111
350 cc	12,7 cm	9,4 cm	6,5 cm	354-6112
450 cc	14,0 cm	10,2 cm	7,1 cm	354-6113
550 cc	15,0 cm	10,9 cm	7,4 cm	354-6114
650 cc	15,7 cm	11,2 cm	7,9 cm	354-6115
750 cc	16,5 cm	11,9 cm	8,1 cm	354-6116

**SILTEX™ CONTOUR PROFILE™**  
 Tissue Expanders, Medium Height,  
 Style 6200\*


Volume	Width	Height	Projection	Catalog #
275 cc	10,7 cm	9,3 cm	6,2 cm	354-6211
350 cc	11,7 cm	10,0 cm	6,6 cm	354-6212
450 cc	12,7 cm	10,8 cm	7,0 cm	354-6213
550 cc	13,5 cm	11,7 cm	7,4 cm	354-6214
650 cc	14,6 cm	12,6 cm	7,6 cm	354-6215
800 cc	15,6 cm	13,3 cm	8,0 cm	354-6216

**SILTEX™ CONTOUR PROFILE™**  
 Tissue Expanders, Tall Height,  
 Style 6300\*


Volume	Width	Height	Projection	Catalog #
250 cc	10,1 cm	10,7 cm	5,6 cm	354-6311
350 cc	11,3 cm	11,8 cm	6,0 cm	354-6312
450 cc	12,3 cm	12,9 cm	6,5 cm	354-6313
550 cc	13,2 cm	13,8 cm	6,9 cm	354-6314
650 cc	14,0 cm	14,6 cm	7,3 cm	354-6315
750 cc	14,6 cm	15,3 cm	7,6 cm	354-6316
850 cc	15,4 cm	15,9 cm	7,9 cm	354-6317

\*Includes CENTERSCOPE™ Magnetic Detection Device


PLEASE NOTE: CONTOUR PROFILE™ Tissue Expanders are not MRI compatible

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

## Other Tissue Expanders

## Smooth Round Tissue Expanders




Device Volume	Diameter	Projection	Catalog #
400 cc	11,3 cm	5,9 cm	350-4305M
550 cc	13,4 cm	5,8 cm	350-4307M
700 cc	14,6 cm	6,7 cm	350-4304M
850 cc	15,0 cm	7,2 cm	350-4308M
1000 cc	15,5 cm	7,8 cm	350-4306M

## Smooth Elliptical Tissue Expanders




Device Volume	Length	Width	Nominal Projection	TEMPORARY OVEREXPANSION		Catalog #
				Projection	Max. Volume	
25 cc	6,0 cm	4,0 cm	2,3 cm	4,5 cm	75 cc	350-5301M
50 cc	8,0 cm	5,0 cm	2,5 cm	5,3 cm	150 cc	350-5302M
75 cc	10,0 cm	6,0 cm	2,1 cm	6,0 cm	225 cc	350-5303M
150 cc	12,0 cm	7,0 cm	3,2 cm	7,2 cm	450 cc	350-5304M
250 cc	14,0 cm	8,0 cm	4,5 cm	9,0 cm	750 cc	350-5305M
125 cc	15,0 cm	5,0 cm	3,3 cm	6,5 cm	375 cc	350-5306M

## Smooth Crescent Tissue Expanders



Device Volume	Length	Width	Projection	Catalog #
50 cc	9,0 cm	5,7 cm	2,7 cm	350-4312M
100 cc	11,0 cm	6,9 cm	3,7 cm	350-4313M
200 cc	14,6 cm	7,5 cm	4,3 cm	350-4314M
300 cc	16,5 cm	10,3 cm	4,7 cm	350-4315M
400 cc	18,2 cm	11,4 cm	5,2 cm	350-4316M
500 cc	19,6 cm	12,4 cm	5,5 cm	350-4317M

## Smooth Rectangle Tissue Expanders



Device Volume	Length	Width	Projection	Catalog #
50 cc	6,9 cm	3,2 cm	3,2 cm	350-4309M
100 cc	6,6 cm	5,1 cm	5,6 cm	350-4301M
250 cc	10,4 cm	7,3 cm	5,6 cm	350-4302M
400 cc	10,6 cm	9,3 cm	6,6 cm	350-4303M
550 cc	13,7 cm	7,1 cm	7,7 cm	350-4310M
700 cc	15,0 cm	8,0 cm	7,9 cm	350-4311M



## Round Saline Breast Implants

SILTEX™ Round Breast Implants,  
Moderate Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer <sup>1</sup>
Volume	Diameter	Projection	Volume	Diameter	Projection		
125 cc	9,5 cm	3,0 cm	150 cc			354-2610	351-125SZ
150 cc	10,0 cm	3,1 cm	175 cc			354-2615	
175 cc	10,6 cm	3,3 cm	200 cc			354-2620	351-175SZ
200 cc	11,0 cm	3,4 cm	225 cc			354-2625	
225 cc	11,5 cm	3,5 cm	250 cc			354-2630	351-225SZ
250 cc	11,9 cm	3,6 cm	275 cc			354-2635	
275 cc	12,3 cm	3,7 cm	300 cc			354-2640	351-275SZ
300 cc	12,6 cm	3,7 cm	325 cc			354-2645	
325 cc	13,0 cm	3,8 cm	375 cc			354-2650	351-325SZ
350 cc	13,3 cm	3,9 cm	400 cc			354-2655	
375 cc	13,6 cm	4,0 cm	425 cc			354-2660	351-375SZ
425 cc	14,2 cm	4,1 cm	475 cc			354-2670	351-425SZ
475 cc	14,8 cm	4,2 cm	525 cc			354-2680	351-475SZ

Smooth Round Breast Implants,  
Moderate Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer <sup>1</sup>
Volume	Diameter	Projection	Volume	Diameter	Projection		
125 cc	9,5 cm	3,0 cm	150 cc	9,1 cm	3,4 cm	350-1610	351-125SZ
150 cc	10,0 cm	3,1 cm	175 cc	9,6 cm	3,7 cm	350-1615	
175 cc	10,6 cm	3,3 cm	200 cc	10,0 cm	3,9 cm	350-1620	351-175SZ
200 cc	11,0 cm	3,4 cm	225 cc	10,4 cm	3,9 cm	350-1625	
225 cc	11,5 cm	3,5 cm	250 cc	10,9 cm	4,0 cm	350-1630	351-225SZ
250 cc	11,9 cm	3,6 cm	275 cc	11,2 cm	4,1 cm	350-1635	
275 cc	12,3 cm	3,7 cm	300 cc	11,7 cm	4,2 cm	350-1640	351-275SZ
300 cc	12,6 cm	3,7 cm	325 cc	12,1 cm	4,3 cm	350-1645	
325 cc	13,0 cm	3,8 cm	375 cc	12,3 cm	4,5 cm	350-1650	351-325SZ
350 cc	13,3 cm	3,9 cm	400 cc	12,6 cm	4,7 cm	350-1655	
375 cc	13,6 cm	4,0 cm	425 cc	12,8 cm	4,8 cm	350-1660	351-375SZ
425 cc	14,2 cm	4,1 cm	475 cc	13,5 cm	4,9 cm	350-1670	351-425SZ
475 cc	14,8 cm	4,2 cm	525 cc	14,2 cm	4,9 cm	350-1680	351-475SZ
525 cc	15,0 cm	4,2 cm	575 cc	14,5 cm	5,1 cm	350-1685	
575 cc	15,0 cm	4,5 cm	625 cc	14,8 cm	5,1 cm	350-1690	
625 cc	15,2 cm	4,6 cm	700 cc	14,9 cm	5,6 cm	350-1695	
700 cc	15,6 cm	4,9 cm	775 cc	14,9 cm	5,8 cm	350-1697	

<sup>1</sup>Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

## Round Saline Breast Implants

Smooth Round Breast Implants,  
Moderate Plus Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer <sup>1</sup>
Volume	Diameter	Projection	Volume	Diameter	Projection		
175 cc	9,5 cm	3,5 cm	210 cc	9,3 cm	4,2 cm	350-2175	351-2175SZ
200 cc	10,0 cm	3,7 cm	240 cc	9,7 cm	4,4 cm	350-2200	351-2200SZ
225 cc	10,4 cm	3,8 cm	270 cc	10,1 cm	4,6 cm	350-2225	351-2225SZ
250 cc	10,8 cm	4,0 cm	300 cc	10,5 cm	4,7 cm	350-2250	351-2250SZ
275 cc	11,0 cm	4,1 cm	330 cc	10,8 cm	4,9 cm	350-2275	351-2275SZ
300 cc	11,5 cm	4,3 cm	360 cc	11,2 cm	5,0 cm	350-2300	351-2300SZ
325 cc	11,9 cm	4,4 cm	390 cc	11,5 cm	5,2 cm	350-2325	351-2325SZ
350 cc	12,1 cm	4,5 cm	420 cc	11,7 cm	5,3 cm	350-2350	351-2350SZ
375 cc	12,3 cm	4,6 cm	450 cc	12,0 cm	5,4 cm	350-2375	351-2375SZ
400 cc	12,6 cm	4,7 cm	480 cc	12,3 cm	5,5 cm	350-2400	351-2400SZ
425 cc	12,9 cm	4,8 cm	510 cc	12,5 cm	5,6 cm	350-2425	351-2425SZ
450 cc	13,0 cm	4,9 cm	540 cc	12,8 cm	5,7 cm	350-2450	351-2450SZ
475 cc	13,3 cm	5,0 cm	570 cc	13,0 cm	5,8 cm	350-2475	351-2475SZ
500 cc	13,6 cm	5,1 cm	600 cc	13,2 cm	5,9 cm	350-2500	351-2500SZ
550 cc	14,0 cm	5,3 cm	660 cc	13,7 cm	6,1 cm	350-2550	351-2550SZ
600 cc	14,5 cm	5,5 cm	720 cc	14,1 cm	6,3 cm	350-2600	351-2600SZ
650 cc	14,8 cm	5,6 cm	780 cc	14,4 cm	6,5 cm	350-2650	351-2650SZ
700 cc	15,2 cm	5,8 cm	840 cc	14,8 cm	6,6 cm	350-2700	351-2700SZ
750 cc	15,6 cm	5,9 cm	900 cc	15,1 cm	6,8 cm	350-2750	351-2750SZ
800 cc	16,1 cm	6,0 cm	960 cc	15,5 cm	7,0 cm	350-2800	351-2800SZ

Smooth Round Breast Implants,  
High Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #	Sterile Saline Sizer <sup>1</sup>
Volume	Diameter	Projection	Volume	Diameter	Projection		
170 cc	9,0 cm	4,0 cm	200 cc	8,8 cm	4,9 cm	350-3170	
190 cc	9,3 cm	4,1 cm	225 cc	9,1 cm	5,1 cm	350-3190	351-3190SZ
210 cc	9,6 cm	4,2 cm	250 cc	9,5 cm	5,3 cm	350-3210	
230 cc	10,0 cm	4,3 cm	275 cc	9,8 cm	5,5 cm	350-3230	351-3230SZ
250 cc	10,2 cm	4,5 cm	300 cc	10,0 cm	5,6 cm	350-3250	
270 cc	10,4 cm	4,6 cm	325 cc	10,2 cm	5,8 cm	350-3270	351-3270SZ
290 cc	10,8 cm	4,7 cm	350 cc	10,5 cm	5,9 cm	350-3290	351-3290SZ
310 cc	11,0 cm	4,8 cm	375 cc	10,7 cm	6,0 cm	350-3310	351-3310SZ
330 cc	11,3 cm	4,8 cm	400 cc	11,0 cm	6,2 cm	350-3330	351-3330SZ
380 cc	11,7 cm	5,2 cm	450 cc	11,4 cm	6,4 cm	350-3380	351-3380SZ
420 cc	12,0 cm	5,4 cm	500 cc	11,7 cm	6,7 cm	350-3420	351-3420SZ
460 cc	12,4 cm	5,5 cm	550 cc	12,1 cm	6,9 cm	350-3460	
500 cc	12,8 cm	5,6 cm	600 cc	12,4 cm	7,1 cm	350-3500	351-3500SZ
560 cc	13,4 cm	5,7 cm	675 cc	13,1 cm	7,3 cm	350-3560	
630 cc	13,8 cm	5,9 cm	750 cc	13,4 cm	7,4 cm	350-3630	351-3630SZ

<sup>1</sup>Single Use Only

## CONTOUR PROFILE™ Saline Breast Implants

SILTEX™ CONTOUR PROFILE™  
Breast Implants, Moderate Profile

Volume	NOMINAL FILL			MAXIMUM FILL	
	Width	Height	Projection	Volume	Catalog #
175 cc	10,2 cm	8,5 cm	4,3 cm	200 cc	354-2910
225 cc	11,2 cm	9,4 cm	4,5 cm	250 cc	354-2911
275 cc	12,2 cm	10,3 cm	4,5 cm	300 cc	354-2912
350 cc	13,1 cm	11,0 cm	4,7 cm	400 cc	354-2913
425 cc	13,8 cm	11,6 cm	5,2 cm	475 cc	354-2914
525 cc	14,9 cm	12,5 cm	5,6 cm	575 cc	354-2915


SILTEX™ CONTOUR PROFILE™  
Breast Implants, High Profile

Volume	NOMINAL FILL			MAXIMUM FILL		Sterile Saline Sizer <sup>1</sup>
	Width	Height	Projection	Volume	Catalog #	
275 cc	11,5 cm	9,5 cm	5,1 cm	300 cc	354-2711	351-0711SZ
350 cc	12,3 cm	10,5 cm	5,3 cm	400 cc	354-2712	351-0712SZ
450 cc	13,2 cm	11,0 cm	6,1 cm	500 cc	354-2713	351-0713SZ
550 cc	14,0 cm	11,9 cm	6,4 cm	600 cc	354-2714	351-0714SZ
650 cc	15,0 cm	12,7 cm	6,6 cm	725 cc	354-2715	351-0715SZ

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above


SPECTRUM™  
Post-Operatively Adjustable Saline Breast Implants

## Smooth Round SPECTRUM™ Implants




Temp. Volume	Min. Volume	NOMINAL FILL			MAXIMUM FILL			Catalog #
		Final Volume	Diameter	Projection	Final Volume	Diameter	Projection	
125 cc	105 cc	125 cc	9,4 cm	2,7 cm	150 cc	9,3 cm	3,2 cm	350-1410
175 cc	150 cc	175 cc	10,3 cm	3,1 cm	210 cc	10,2 cm	3,8 cm	350-1420
225 cc	190 cc	225 cc	10,8 cm	3,5 cm	270 cc	10,9 cm	4,3 cm	350-1430
275 cc	230 cc	275 cc	11,6 cm	3,8 cm	330 cc	11,5 cm	4,6 cm	350-1440
325 cc	275 cc	325 cc	12,4 cm	4,0 cm	390 cc	12,1 cm	4,9 cm	350-1450
375 cc	320 cc	375 cc	13,5 cm	4,1 cm	450 cc	13,0 cm	4,9 cm	350-1460
425 cc	360 cc	425 cc	14,0 cm	4,2 cm	510 cc	13,7 cm	4,8 cm	350-1470
475 cc	405 cc	475 cc	14,3 cm	4,4 cm	570 cc	13,8 cm	5,4 cm	350-1480
525 cc	450 cc	525 cc	15,0 cm	4,4 cm	630 cc	14,4 cm	5,6 cm	350-1485
575 cc	490 cc	575 cc	14,7 cm	4,8 cm	690 cc	14,2 cm	6,0 cm	350-1490

## SILTEX™ Round SPECTRUM™ Implants



Temp. Volume	Min. Volume	NOMINAL FILL			MAXIMUM FILL			Catalog #
		Final Volume	Diameter	Projection	Final Volume	Diameter	Projection	
125 cc	105 cc	125 cc	9,4 cm	3,0 cm	150 cc	9,2 cm	3,5 cm	354-2410M
175 cc	150 cc	175 cc	10,5 cm	3,1 cm	210 cc	10,4 cm	3,9 cm	354-2420M
225 cc	190 cc	225 cc	11,3 cm	3,3 cm	270 cc	11,1 cm	4,2 cm	354-2430M
275 cc	235 cc	275 cc	12,0 cm	3,6 cm	330 cc	11,4 cm	4,4 cm	354-2440M
325 cc	275 cc	325 cc	12,9 cm	3,7 cm	390 cc	12,6 cm	4,5 cm	354-2450M
375 cc	320 cc	375 cc	13,4 cm	3,9 cm	450 cc	13,1 cm	4,8 cm	354-2460M
425 cc	360 cc	425 cc	14,0 cm	4,2 cm	510 cc	13,6 cm	5,0 cm	354-2470M
475 cc	400 cc	475 cc	14,5 cm	4,5 cm	570 cc	14,2 cm	5,7 cm	354-2480M

SILTEX™ CONTOUR PROFILE™  
Expanders/Implants

Temp. Volume	Min. Volume	NOMINAL FILL				MAXIMUM FILL				Catalog #
		Final Volume	Width	Height	Proj.	Final Volume	Width	Height	Proj.	
275 cc	235 cc	275 cc	11,5 cm	9,5 cm	5,1 cm	330 cc	11,1 cm	9,4 cm	6,3 cm	354-2511
350 cc	300 cc	350 cc	12,3 cm	10,5 cm	5,3 cm	420 cc	11,9 cm	10,1 cm	6,9 cm	354-2512
450 cc	380 cc	450 cc	13,2 cm	11,0 cm	6,1 cm	540 cc	12,7 cm	11,0 cm	7,5 cm	354-2513
550 cc	470 cc	550 cc	14,0 cm	11,9 cm	6,4 cm	660 cc	13,7 cm	11,9 cm	7,9 cm	354-2514
650 cc	550 cc	650 cc	15,0 cm	12,7 cm	6,6 cm	780 cc	14,5 cm	12,5 cm	8,2 cm	354-2515

Sterile Saline Breast Implant Sizers<sup>†</sup>Round Saline Breast Implant Sizers,  
Moderate Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
125 cc	9,5 cm	3,0 cm	150 cc	9,1 cm	3,4 cm	351-125SZ
175 cc	10,6 cm	3,3 cm	200 cc	10,0 cm	3,9 cm	351-175SZ
225 cc	11,5 cm	3,5 cm	250 cc	10,9 cm	4,0 cm	351-225SZ
275 cc	12,3 cm	3,7 cm	300 cc	11,7 cm	4,2 cm	351-275SZ
325 cc	13,0 cm	3,8 cm	375 cc	12,3 cm	4,5 cm	351-325SZ
375 cc	13,6 cm	4,0 cm	425 cc	12,8 cm	4,8 cm	351-375SZ
425 cc	14,2 cm	4,1 cm	475 cc	13,5 cm	4,9 cm	351-425SZ
475 cc	14,8 cm	4,2 cm	525 cc	14,2 cm	4,9 cm	351-475SZ

Round Saline Breast Implant Sizers,  
Moderate Plus Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
175 cc	9,5 cm	3,5 cm	210 cc	9,3 cm	4,2 cm	351-2175SZ
200 cc	10,0 cm	3,7 cm	240 cc	9,7 cm	4,4 cm	351-2200SZ
225 cc	10,4 cm	3,8 cm	270 cc	10,1 cm	4,6 cm	351-2225SZ
250 cc	10,8 cm	4,0 cm	300 cc	10,5 cm	4,7 cm	351-2250SZ
275 cc	11,0 cm	4,1 cm	330 cc	10,8 cm	4,9 cm	351-2275SZ
300 cc	11,5 cm	4,3 cm	360 cc	11,2 cm	5,0 cm	351-2300SZ
325 cc	11,9 cm	4,4 cm	390 cc	11,5 cm	5,2 cm	351-2325SZ
350 cc	12,1 cm	4,5 cm	420 cc	11,7 cm	5,3 cm	351-2350SZ
375 cc	12,3 cm	4,6 cm	450 cc	12,0 cm	5,4 cm	351-2375SZ
400 cc	12,6 cm	4,7 cm	480 cc	12,3 cm	5,5 cm	351-2400SZ
425 cc	12,9 cm	4,8 cm	510 cc	12,5 cm	5,6 cm	351-2425SZ
450 cc	13,0 cm	4,9 cm	540 cc	12,8 cm	5,7 cm	351-2450SZ
475 cc	13,3 cm	5,0 cm	570 cc	13,0 cm	5,8 cm	351-2475SZ
500 cc	13,6 cm	5,1 cm	600 cc	13,2 cm	5,9 cm	351-2500SZ
550 cc	14,0 cm	5,3 cm	660 cc	13,7 cm	6,1 cm	351-2550SZ
600 cc	14,5 cm	5,5 cm	720 cc	14,1 cm	6,3 cm	351-2600SZ
650 cc	14,8 cm	5,6 cm	780 cc	14,4 cm	6,5 cm	351-2650SZ
700 cc	15,2 cm	5,8 cm	840 cc	14,8 cm	6,6 cm	351-2700SZ
750 cc	15,6 cm	5,9 cm	900 cc	15,1 cm	6,8 cm	351-2750SZ
800 cc	16,1 cm	6,0 cm	960 cc	15,5 cm	7,0 cm	351-2800SZ

Sterile Saline Implant Sizers<sup>†</sup> CONT'DRound Saline Breast Implant Sizers,  
High Profile

NOMINAL FILL			MAXIMUM FILL			Catalog #
Volume	Diameter	Projection	Volume	Diameter	Projection	
190 cc	9,3 cm	4,1 cm	225 cc	9,1 cm	5,1 cm	351-3190SZ
230 cc	10,0 cm	4,3 cm	275 cc	9,8 cm	5,5 cm	351-3230SZ
270 cc	10,4 cm	4,6 cm	325 cc	10,2 cm	5,8 cm	351-3270SZ
290 cc	10,8 cm	4,7 cm	350 cc	10,5 cm	5,9 cm	351-3290SZ
310 cc	11,0 cm	4,8 cm	375 cc	10,7 cm	6,0 cm	351-3310SZ
330 cc	11,3 cm	4,8 cm	400 cc	11,0 cm	6,2 cm	351-3330SZ
380 cc	11,7 cm	5,2 cm	450 cc	11,4 cm	6,4 cm	351-3380SZ
420 cc	12,0 cm	5,4 cm	500 cc	11,7 cm	6,7 cm	351-3420SZ
500 cc	12,8 cm	5,6 cm	600 cc	12,4 cm	7,1 cm	351-3500SZ
630 cc	13,8 cm	5,9 cm	750 cc	13,4 cm	7,4 cm	351-3630SZ

CONTOUR PROFILE™ Saline  
Breast Implant Sizers, High Profile

Volume	Width	Height	Projection	Catalog #
275 + 25 cc	11,5 cm	9,5 cm	5,1 cm	351-0711SZ
350 + 50 cc	12,3 cm	10,5 cm	5,3 cm	351-0712SZ
450 + 50 cc	13,2 cm	11,0 cm	6,1 cm	351-0713SZ
550 + 50 cc	14,0 cm	11,9 cm	6,4 cm	351-0714SZ
650 + 75 cc	15,0 cm	12,7 cm	6,6 cm	351-0715SZ

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above


Accessories

Aseptic Transfer System

Description	Catalog #
*Aseptic Transfer System	350-8400

Carton contains 4 sterile kits, each kit includes:

- 60 cc syringe (sterile, disposable)
- Piercing spike
- 48" tubing (saline container to syringe)
- Dual check valve
- 10" tubing (syringe to implant fill tube)
- Luer connector (attaches to implant fill tube)



Injection Port Detector

Description	Catalog #
CENTERSCOPE™ Magnetic Detection Device	350-4402

Injection Dome

Description	Catalog #
Dome Pack Set	350-DOMPK
Can be used with SPECTRA™ Adjustable Gel Implant.	

Important Safety Information

MemoryGel™ Breast Implants and MENTOR® Saline-Filled Breast Implants are indicated for breast augmentation in women (at least 18 years old) or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body with existing cancer or precancer of their breast who have not received adequate treatment for those conditions or are pregnant or nursing.

Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery. The most common complications with the MemoryGel™ Breast Implants include reoperation, capsular contracture, asymmetry, and breast pain. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.

The most common complications with MENTOR® Saline-Filled Breast Implants include reoperation, implant removal, capsular contracture, wrinkling, breast pain and deflation.

Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.

CONTOUR PROFILE™ Tissue Expanders and Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy. These expanders are intended for temporary subcutaneous or submuscular implantation and are not intended for use beyond six months. Do not use the CONTOUR PROFILE™ Tissue Expander and Smooth Tissue Expander with Remote Dome in patients where an MRI may be needed. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas.

For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, which include MENTOR® Saline-Filled Breast Implants, MemoryGel™ Breast Implants and Tissue Expanders, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at [www.mentorwillc.com](http://www.mentorwillc.com).

\*Not available for sale in the European Union

## Corporate Headquarters

Santa Barbara, CA 93111 USA

Tel: +1 805 879 6000

Fax: +1 805 967 7108

www.mentorwllc.com

## Operations Europe, Middle East, Africa

### Germany, Austria & Switzerland

Mentor Deutschland GmbH

Tel: +49 811 600 50 0

Fax: +49 811 600 50 11

### Italy, Greece & Israel

Johnson & Johnson Medical S.p.a.

Tel: +39 06 911 94 500

Fax: +39 02 644 10 563

### France & Benelux

ETHICON France, S.A.S.

Tel: +33 1 55 00 22 45

Fax: +33 1 55 00 28 11

### Croatia, Macedonia, Serbia, Slovenia & Bulgaria

Johnson & Johnson S.E. d.o.o.

Tel: +386 1 40 11 800

Fax: +386 1 40 11 801

### Poland, Lithuania, Latvia & Estonia

Johnson & Johnson Poland sp. z o.o.

Tel: +48 22 237 83 05

Fax: +48 22 237 83 53

### South Africa

Johnson & Johnson Medical (Pty) Ltd

Tel: +27 11 265 1135

Fax: +27 11 265 1114

### Middle East, Egypt & Pakistan

Johnson & Johnson Middle East Inc

Tel: +971 4 429 72 00

Fax: +971 4 299 43 30

### Spain & Portugal

Mentor Medical Systems Iberica, S.L

Tel: +34 91 562 2700

Fax: +34 91 728 3367

### United Kingdom & Ireland

Johnson & Johnson Medical Ltd

Tel: +44 (0)800 973 082

Fax: +44 (0)134 486 4122

### Russia & CIS

Johnson & Johnson LLC

Tel: +7 495 580 7777

Fax: +7 495 580 7878

### Czech Republic, Slovak Republic, Hungary & Romania

Johnson & Johnson s.r.o.

Tel: +420 227 012 371

Fax: +420 227 012 300

### Sweden, Denmark, Norway, Finland & Iceland

Johnson & Johnson AB

Tel: +46 8 626 2200

Fax: +46 8 754 5850

### Turkey

Johnson & Johnson Medical Turkey

Tel: +90 216 538 21 21

Fax: +90 216 538 21 30

### Maghreb & Sub Saharan Africa

Johnson & Johnson Professional Export

Tel: +44 (0)1506 594 500

Fax: +44 (0)1506 594 731



# IMPLANT CATALOGUE

EFFECTIVE JANUARY 2012



**PERTHESE™**  
BREAST IMPLANTS



Make life more beautiful™

NOT FOR DISTRIBUTION IN USA  
APPROVED FOR EMEA REGION

© Mentor Worldwide LLC 2012 0907007 Rev C

Based on the high level of manufacturing standards and quality testing, Perouse Plastique was the first French implant manufacturer given the authorization to sell silicone gel products into the local market after a ban period of several years. The company launched its trademark **PERTHESE™ Breast Implant** in 1990.

During the following years, Perouse Plastique was not only able to continuously strengthen its market position in France but also in many other EU countries as well as in Latin America. With the acquisition through Mentor in 2007, more surgeons worldwide now have the opportunity to access this unique range of products adapted to the needs of medical professionals AND their patients in both Aesthetic and Reconstructive procedures.

The **PERTHESE™ Breast Implant** is a sterile gel-filled implant made of a trilaminar silicone envelope. The silicone envelope consists of an internal and external layer made with a high mechanical resistance medical grade silicone elastomer and an intermediate barrier layer to significantly reduce gel bleed.

The **PERTHESE™ Implant** range is available in two lines of breast implants:

### **PERTHESE™ Breast Implants**

Round breast implants filled with a soft and cohesive silicone gel, available either with a micro-textured surface or with a smooth surface.

Product Profile	Micro-textured	Smooth
TX, Low Profile	X	
MX, Moderate Profile	X	
MP, Moderate Plus Profile	X	
540, High Profile	X	
550, Ultra High Profile	X	
700T, Low Profile		X
510T, High Profile		X

### **PERTHESE™ Esthea Breast Implants**

Anatomical breast implants filled with a cohesive silicone gel with a micro-textured surface, available in different base shapes.

Product Profile	Micro-textured	Smooth
ELP-EHP-ESHP	X	
AX	X	
EFP	X	
EMP-EMHP	X	

All **PERTHESE™ Silicone Gel-filled Breast Implants** are covered by **PERTHESE™ Quality Assurance Program warranty** or **PatientCare Promise limited warranty\***

### **PERTHESE™ Silicone Sheetings**

**PERTHESE™ Silicone Sheeting** is a flexible, translucent silicone elastomer sheeting material designed for medical and laboratory applications. It is made from an enhanced tear resistant elastomer that consists of dimethyl and methyl vinyl siloxane copolymers.

† In case of eligibility, depending on the date of surgery, the implant will fall under one of the above warranties.

## CONTENTS

### **PERTHESE™ Breast Implants Cohesive Silicone Gel Breast Implants, Round Shapes. .... 1**

PERTHESE™ Breast Implants TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implants MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implants MP, Moderate Plus Profile, Micro-Textured  
PERTHESE™ Breast Implants 540, High Profile, Micro-Textured  
PERTHESE™ Breast Implants 550, Ultra High Profile, Micro-Textured  
PERTHESE™ Breast Implants 700T, Low Profile, Smooth  
PERTHESE™ Breast Implants 510T, High Profile, Smooth  
PERTHESE™ Breast Implants Summary Chart

### **PERTHESE™ Esthea Breast Implants Cohesive Silicone Gel Breast Implants, Anatomical Shapes. .... 7**

PERTHESE™ Esthea Breast Implants ELP - EHP - ESHP, Micro-Textured  
PERTHESE™ Esthea Breast Implants AX, Micro-Textured  
PERTHESE™ Esthea Breast Implants EFP, Micro-Textured  
PERTHESE™ Esthea Breast Implants EMP - EMHP, Micro-Textured  
MENTOR® CPX2™ Tissue Expanders / PERTHESE™ Esthea Breast  
Implants Sizing Guide – Two-Stage Breast Reconstruction

### **PERTHESE™ Sterile Silicone Breast Implant Sizers, Round Shapes ..... 11**

PERTHESE™ Breast Implant Sizers TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers 540, High Profile, Micro-Textured

### **PERTHESE™ Sterile Inflatable Breast Implant Sizers, Round Shapes ..... 13**

PERTHESE™ Breast Implant Sizers TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers 540, High Profile, Micro-Textured

### **PERTHESE™ Esthea Sterile Silicone Breast Implant Sizers, Anatomical Shapes. .... 15**

PERTHESE™ Breast Esthea Implant Sizers ELP - EHP - ESHP,  
Micro-Textured  
PERTHESE™ Breast Esthea Implant Sizers AX, Micro-Textured  
PERTHESE™ Breast Esthea Implant Sizers EMP, Micro-Textured

### **PERTHESE™ Esthea Sterile Inflatable Breast Implant Sizers, Anatomical Shapes. .... 17**

PERTHESE™ Breast Implant Sizers ELP - EHP - ESHP, Micro-Textured  
PERTHESE™ Breast Implant Sizers AX, Micro-Textured  
PERTHESE™ Breast Implant Sizers EMP, Micro-Textured

### **PERTHESE™ Silicone Sheetings..... 19**

### **Contact Information ..... 20**



Cohesive Silicone Gel Breast Implants,  
Round Shapes

PERTHESE™ Breast Implants TX,  
Low Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
100 cc	9,1 cm	2,2 cm	<b>TX 100</b>
125 cc	9,7 cm	2,4 cm	<b>TX 125</b>
150 cc	10,2 cm	2,6 cm	<b>TX 150</b>
170 cc	10,8 cm	2,7 cm	<b>TX 170</b>
185 cc	11,1 cm	2,8 cm	<b>TX 185</b>
200 cc	11,4 cm	2,9 cm	<b>TX 200</b>
220 cc	11,8 cm	2,9 cm	<b>TX 220</b>
240 cc	12,2 cm	3,0 cm	<b>TX 240</b>
260 cc	12,3 cm	3,1 cm	<b>TX 260</b>
280 cc	12,5 cm	3,2 cm	<b>TX 280</b>
300 cc	12,8 cm	3,2 cm	<b>TX 300</b>
325 cc	13,1 cm	3,3 cm	<b>TX 325</b>
350 cc	13,5 cm	3,4 cm	<b>TX 350</b>

PERTHESE™ Breast Implants MX,  
Moderate Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
100 cc	8,7 cm	2,8 cm	<b>MX 100</b>
125 cc	9,2 cm	3,0 cm	<b>MX 125</b>
150 cc	9,9 cm	3,2 cm	<b>MX 150</b>
175 cc	10,4 cm	3,3 cm	<b>MX 175</b>
200 cc	11,0 cm	3,4 cm	<b>MX 200</b>
225 cc	11,5 cm	3,5 cm	<b>MX 225</b>
250 cc	12,0 cm	3,6 cm	<b>MX 250</b>
275 cc	12,2 cm	3,7 cm	<b>MX 275</b>
300 cc	12,5 cm	3,8 cm	<b>MX 300</b>
325 cc	12,8 cm	4,0 cm	<b>MX 325</b>
350 cc	13,0 cm	4,2 cm	<b>MX 350</b>
375 cc	13,2 cm	4,4 cm	<b>MX 375</b>
400 cc	13,5 cm	4,5 cm	<b>MX 400</b>
450 cc	13,9 cm	4,6 cm	<b>MX 450</b>
500 cc	14,2 cm	5,0 cm	<b>MX 500</b>
550 cc	14,5 cm	5,2 cm	<b>MX 550</b>
600 cc	14,9 cm	5,4 cm	<b>MX 600</b>
650 cc	15,2 cm	5,5 cm	<b>MX 650</b>
700 cc	15,6 cm	5,6 cm	<b>MX 700</b>
800 cc	16,1 cm	6,0 cm	<b>MX 800</b>

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Cohesive Silicone Gel Breast Implants,  
Round Shapes CONT'D

PERTHESE™ Breast Implants MP,  
Moderate Plus Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
115 cc	8,7 cm	3,1 cm	<b>MP 115</b>
135 cc	9,2 cm	3,2 cm	<b>MP 135</b>
165 cc	9,9 cm	3,4 cm	<b>MP 165</b>
190 cc	10,5 cm	3,6 cm	<b>MP 190</b>
225 cc	11,0 cm	3,9 cm	<b>MP 225</b>
260 cc	11,6 cm	4,1 cm	<b>MP 260</b>
305 cc	12,1 cm	4,4 cm	<b>MP 305</b>
320 cc	12,3 cm	4,5 cm	<b>MP 320</b>
345 cc	12,5 cm	4,6 cm	<b>MP 345</b>
380 cc	12,7 cm	4,9 cm	<b>MP 380</b>
410 cc	13,0 cm	5,0 cm	<b>MP 410</b>
440 cc	13,2 cm	5,3 cm	<b>MP 440</b>
475 cc	13,5 cm	5,4 cm	<b>MP 475</b>
510 cc	13,8 cm	5,5 cm	<b>MP 510</b>
580 cc	14,3 cm	5,6 cm	<b>MP 580</b>
635 cc	14,4 cm	5,9 cm	<b>MP 635</b>
685 cc	14,9 cm	6,1 cm	<b>MP 685</b>
735 cc	15,2 cm	6,3 cm	<b>MP 735</b>
800 cc	15,7 cm	6,4 cm	<b>MP 800</b>
910 cc	16,2 cm	6,8 cm	<b>MP 910</b>

PERTHESE™ Breast Implants 540,  
High Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
75 cc	8,1 cm	2,3 cm	<b>540-075</b>
100 cc	8,2 cm	3,0 cm	<b>540-100</b>
125 cc	8,8 cm	3,2 cm	<b>540-125</b>
150 cc	9,2 cm	3,4 cm	<b>540-150</b>
175 cc	9,9 cm	3,6 cm	<b>540-175</b>
200 cc	10,1 cm	3,8 cm	<b>540-200</b>
230 cc	10,4 cm	4,0 cm	<b>540-230</b>
260 cc	11,0 cm	4,3 cm	<b>540-260</b>
280 cc	11,1 cm	4,4 cm	<b>540-280</b>
300 cc	11,2 cm	4,5 cm	<b>540-300</b>
325 cc	11,7 cm	4,6 cm	<b>540-325</b>
350 cc	11,9 cm	4,8 cm	<b>540-350</b>
380 cc	12,0 cm	4,9 cm	<b>540-380</b>
400 cc	12,1 cm	4,9 cm	<b>540-400</b>
440 cc	12,7 cm	5,3 cm	<b>540-440</b>
500 cc	13,2 cm	5,4 cm	<b>540-500</b>
550 cc	13,25 cm	5,9 cm	<b>540-550</b>
600 cc	13,8 cm	6,0 cm	<b>540-600</b>
650 cc	14,2 cm	6,1 cm	<b>540-650</b>
850 cc	15,5 cm	6,7 cm	<b>540-850</b>

Cohesive Silicone Gel Breast Implants,  
Round Shapes CONT'D

**PERTHESE™ Breast Implants 550,  
Ultra High Profile, Micro-Textured**



Volume	Diameter	Projection	Catalog #
95 cc	7,5 cm	3,5 cm	<b>550-095</b>
125 cc	8,2 cm	3,8 cm	<b>550-125</b>
150 cc	8,8 cm	3,9 cm	<b>550-150</b>
160 cc	9,0 cm	4,0 cm	<b>550-160</b>
175 cc	9,2 cm	4,1 cm	<b>550-175</b>
195 cc	9,5 cm	4,2 cm	<b>550-195</b>
215 cc	9,7 cm	4,3 cm	<b>550-215</b>
235 cc	10,1 cm	4,5 cm	<b>550-235</b>
265 cc	10,4 cm	4,7 cm	<b>550-265</b>
285 cc	10,8 cm	4,7 cm	<b>550-285</b>
315 cc	11,1 cm	4,9 cm	<b>550-315</b>
330 cc	11,2 cm	5,1 cm	<b>550-330</b>
345 cc	11,3 cm	5,2 cm	<b>550-345</b>
390 cc	11,7 cm	5,4 cm	<b>550-390</b>
410 cc	12,0 cm	5,4 cm	<b>550-410</b>
425 cc	12,1 cm	5,5 cm	<b>550-425</b>
445 cc	12,2 cm	5,7 cm	<b>550-445</b>
485 cc	12,6 cm	6,0 cm	<b>550-485</b>
530 cc	12,9 cm	6,0 cm	<b>550-530</b>
580 cc	13,4 cm	6,1 cm	<b>550-580</b>
630 cc	13,6 cm	6,3 cm	<b>550-630</b>
705 cc	14,1 cm	6,7 cm	<b>550-705</b>

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Cohesive Silicone Gel Breast Implants,  
Round Shapes CONT'D

**PERTHESE™ Breast Implants 700T,  
Low Profile, Smooth**



Volume	Diameter	Projection	Catalog #
120 cc	9,7 cm	2,4 cm	<b>700 T-120</b>
140 cc	10,3 cm	2,4 cm	<b>700 T-140</b>
160 cc	10,8 cm	2,5 cm	<b>700 T-160</b>
180 cc	11,2 cm	2,6 cm	<b>700 T-180</b>
200 cc	11,2 cm	2,8 cm	<b>700 T-200</b>
220 cc	11,9 cm	2,8 cm	<b>700 T-220</b>
240 cc	12,1 cm	2,9 cm	<b>700 T-240</b>
260 cc	12,2 cm	3,1 cm	<b>700 T-260</b>
280 cc	12,6 cm	3,2 cm	<b>700 T-280</b>
300 cc	12,8 cm	3,2 cm	<b>700 T-300</b>
325 cc	13,1 cm	3,3 cm	<b>700 T-325</b>
350 cc	13,5 cm	3,3 cm	<b>700 T-350</b>


**PERTHESE™ Breast Implants 510T,  
High Profile, Smooth**



Volume	Diameter	Projection	Catalog #
80 cc	8,1 cm	2,5 cm	<b>510 T-080</b>
100 cc	9,1 cm	2,7 cm	<b>510 T-100</b>
120 cc	9,3 cm	2,9 cm	<b>510 T-120</b>
140 cc	9,6 cm	3,3 cm	<b>510 T-140</b>
160 cc	9,6 cm	3,4 cm	<b>510 T-160</b>
180 cc	10,0 cm	3,6 cm	<b>510 T-180</b>
200 cc	10,2 cm	3,7 cm	<b>510 T-200</b>
220 cc	10,7 cm	3,9 cm	<b>510 T-220</b>
240 cc	10,7 cm	4,2 cm	<b>510 T-240</b>
270 cc	11,0 cm	4,3 cm	<b>510 T-270</b>
300 cc	11,4 cm	4,4 cm	<b>510 T-300</b>
340 cc	11,7 cm	4,6 cm	<b>510 T-340</b>
380 cc	11,9 cm	4,8 cm	<b>510 T-380</b>
400 cc	12,2 cm	5,0 cm	<b>510 T-400</b>
450 cc	12,6 cm	5,3 cm	<b>510 T-450</b>
500 cc	13,3 cm	5,5 cm	<b>510 T-500</b>



Cohesive Silicone Gel Breast Implants,  
Anatomical Shapes



PERTHESE™ Esthea Breast Implants  
ELP - EHP - ESHP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	<b>ELP 190</b>
225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	<b>EHP 225</b>
235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	<b>ESHP 235</b>
240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	<b>ELP 240</b>
275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	<b>EHP 275</b>
295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	<b>ESHP 295</b>
270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	<b>ELP 270</b>
325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	<b>EHP 325</b>
355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	<b>ESHP 355</b>
310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	<b>ELP 310</b>
375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	<b>EHP 375</b>
415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	<b>ESHP 415</b>
360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	<b>ELP 360</b>
425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	<b>EHP 425</b>
475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	<b>ESHP 475</b>

PERTHESE™ Esthea Breast Implants  
AX, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
190 cc	10,6 cm	10,8 cm	4,3 cm	3,7 cm	5,8 cm	7,8 cm	<b>AX 190</b>
210 cc	10,9 cm	11,2 cm	4,4 cm	3,8 cm	6,0 cm	8,2 cm	<b>AX 210</b>
240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	<b>AX 240</b>
260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	<b>AX 260</b>
285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	<b>AX 285</b>
300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	<b>AX 300</b>
340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	<b>AX 340</b>
380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	<b>AX 380</b>
400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	<b>AX 400</b>
450 cc	13,9 cm	14,7 cm	5,2 cm	4,8 cm	7,5 cm	10,8 cm	<b>AX 450</b>
480 cc	14,5 cm	15,2 cm	5,3 cm	4,9 cm	7,7 cm	11,2 cm	<b>AX 480</b>
500 cc	14,9 cm	15,6 cm	5,3 cm	5,0 cm	7,9 cm	11,5 cm	<b>AX 500</b>

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Cohesive Silicone Gel Breast Implants,  
Anatomical Shapes CONT'D

PERTHESE™ Esthea Breast Implants  
EFP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
95 cc	8,1 cm	7,5 cm	3,4 cm	2,7 cm	4,2 cm	5,4 cm	<b>EFP 095</b>
115 cc	8,6 cm	8,1 cm	3,6 cm	3,0 cm	4,6 cm	5,8 cm	<b>EFP 115</b>
140 cc	9,2 cm	8,5 cm	3,8 cm	3,2 cm	5,0 cm	6,1 cm	<b>EFP 140</b>
180 cc	10,2 cm	9,2 cm	4,1 cm	3,5 cm	5,4 cm	6,6 cm	<b>EFP 180</b>
210 cc	10,7 cm	9,7 cm	4,3 cm	3,7 cm	5,7 cm	6,9 cm	<b>EFP 210</b>
250 cc	11,2 cm	10,2 cm	4,5 cm	3,9 cm	6,1 cm	7,3 cm	<b>EFP 250</b>
285 cc	11,7 cm	10,7 cm	4,7 cm	4,1 cm	6,4 cm	7,7 cm	<b>EFP 285</b>
325 cc	12,2 cm	11,3 cm	4,9 cm	4,2 cm	6,6 cm	8,1 cm	<b>EFP 325</b>
370 cc	12,7 cm	11,8 cm	5,1 cm	4,4 cm	6,9 cm	8,4 cm	<b>EFP 370</b>
415 cc	13,1 cm	12,2 cm	5,3 cm	4,7 cm	7,3 cm	8,7 cm	<b>EFP 415</b>
470 cc	13,6 cm	12,6 cm	5,5 cm	4,9 cm	7,7 cm	9,0 cm	<b>EFP 470</b>
520 cc	14,0 cm	13,1 cm	5,7 cm	5,1 cm	8,0 cm	9,4 cm	<b>EFP 520</b>
575 cc	14,6 cm	13,5 cm	5,9 cm	5,3 cm	8,3 cm	9,6 cm	<b>EFP 575</b>
640 cc	15,1 cm	14,0 cm	6,1 cm	5,6 cm	8,8 cm	9,9 cm	<b>EFP 640</b>
710 cc	15,5 cm	14,2 cm	6,3 cm	5,8 cm	9,1 cm	10,1 cm	<b>EFP 710</b>

PERTHESE™ Esthea Breast Implants  
EMP - EMHP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
80 cc	8,1 cm	6,8 cm	3,2 cm	2,4 cm	3,7 cm	4,9 cm	<b>EMP 080</b>
90 cc	8,1 cm	6,8 cm	3,5 cm	2,4 cm	4,1 cm	5,2 cm	<b>EMHP 090</b>
100 cc	8,6 cm	7,3 cm	3,5 cm	2,7 cm	4,2 cm	5,3 cm	<b>EMP 100</b>
110 cc	8,6 cm	7,3 cm	3,8 cm	2,7 cm	4,6 cm	5,5 cm	<b>EMHP 110</b>
120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	<b>EMP 120</b>
130 cc	9,2 cm	7,7 cm	4,0 cm	2,9 cm	4,9 cm	5,8 cm	<b>EMHP 130</b>
145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	<b>EMP 145</b>
170 cc	10,2 cm	8,5 cm	4,3 cm	3,1 cm	5,3 cm	6,4 cm	<b>EMHP 170</b>
165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	<b>EMP 165</b>
195 cc	10,7 cm	9,1 cm	4,5 cm	3,3 cm	5,6 cm	6,8 cm	<b>EMHP 195</b>
200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	<b>EMP 200</b>
235 cc	11,2 cm	9,6 cm	4,7 cm	3,5 cm	5,9 cm	7,2 cm	<b>EMHP 235</b>
230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	<b>EMP 230</b>
270 cc	11,7 cm	10,1 cm	4,9 cm	3,7 cm	6,3 cm	7,6 cm	<b>EMHP 270</b>
265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	<b>EMP 265</b>
320 cc	12,2 cm	10,6 cm	5,1 cm	3,9 cm	6,6 cm	7,9 cm	<b>EMHP 320</b>
305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	<b>EMP 305</b>
360 cc	12,7 cm	11,1 cm	5,3 cm	4,1 cm	6,9 cm	8,3 cm	<b>EMHP 360</b>
345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	<b>EMP 345</b>
405 cc	13,1 cm	11,6 cm	5,5 cm	4,3 cm	7,2 cm	8,6 cm	<b>EMHP 405</b>
390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	<b>EMP 390</b>
455 cc	13,6 cm	12,0 cm	5,7 cm	4,5 cm	7,6 cm	9,0 cm	<b>EMHP 455</b>
435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	<b>EMP 435</b>
490 cc	14,0 cm	12,5 cm	5,8 cm	4,3 cm	7,7 cm	9,6 cm	<b>EMHP 490</b>
475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	<b>EMP 475</b>
545 cc	14,6 cm	13,0 cm	6,0 cm	4,6 cm	8,0 cm	9,7 cm	<b>EMHP 545</b>
535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	<b>EMP 535</b>
625 cc	15,1 cm	13,3 cm	6,2 cm	5,0 cm	8,6 cm	9,9 cm	<b>EMHP 625</b>
585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	<b>EMP 585</b>
680 cc	15,5 cm	13,6 cm	6,4 cm	5,2 cm	8,9 cm	10,2 cm	<b>EMHP 680</b>

MENTOR® CPX2™ Tissue Expanders /  
PERTHESE™ Esthea Breast Implants Sizing Guide –

Two-Stage Breast Reconstruction

**STAGE 1**

**Selection of the MENTOR® CPX2™ Expander\***

**Style 6100 - Low Height**

Vol.	Width	Height	Proj.	Catalog #
250 cc	11,4 cm	8,1 cm	6,1 cm	354-6111
350 cc	12,7 cm	9,4 cm	6,5 cm	354-6112
450 cc	14,0 cm	10,2 cm	7,1 cm	354-6113
550 cc	15,0 cm	10,9 cm	7,4 cm	354-6114
650 cc	15,7 cm	11,2 cm	7,9 cm	354-6115
750 cc	16,5 cm	11,9 cm	8,1 cm	354-6116

**Style 6200 – Medium Height**

Vol.	Width	Height	Proj.	Catalog #
275 cc	10,7 cm	9,3 cm	6,2 cm	354-6211
350 cc	11,7 cm	10,0 cm	6,6 cm	354-6212
450 cc	12,7 cm	10,8 cm	7,0 cm	354-6213
550 cc	13,5 cm	11,7 cm	7,4 cm	354-6214
650 cc	14,6 cm	12,6 cm	7,6 cm	354-6215
800 cc	15,6 cm	13,3 cm	8,0 cm	354-6216

**Style 6300 - Tall Height**

Vol.	Width	Height	Proj.	Catalog #
250 cc	10,1 cm	10,7 cm	5,6 cm	354-6311
350 cc	11,3 cm	11,8 cm	6,0 cm	354-6312
450 cc	12,3 cm	12,9 cm	6,5 cm	354-6313
550 cc	13,2 cm	13,8 cm	6,9 cm	354-6314
650 cc	14,0 cm	14,6 cm	7,3 cm	354-6315
750 cc	14,6 cm	15,3 cm	7,6 cm	354-6316
850 cc	15,4 cm	15,9 cm	7,9 cm	354-6317

**PERTHESE™ Esthea Breast Implants EMP**

Vol.	Width	Height	Proj.	Catalog #
165 cc	10,7 cm	9,1 cm	4,0 cm	EMP165
200 cc	11,2 cm	9,6 cm	4,3 cm	EMP200
230 cc	11,7 cm	10,1 cm	4,5 cm	EMP230
265 cc	12,2 cm	10,6 cm	4,6 cm	EMP265
305 cc	12,7 cm	11,1 cm	4,9 cm	EMP305
345 cc	13,1 cm	11,6 cm	5,0 cm	EMP345
390 cc	13,6 cm	12,0 cm	5,2 cm	EMP390
435 cc	14,0 cm	12,5 cm	5,5 cm	EMP435
475 cc	14,6 cm	13,0 cm	5,7 cm	EMP475
535 cc	15,1 cm	13,3 cm	5,9 cm	EMP535
585 cc	15,5 cm	13,6 cm	6,1 cm	EMP585

**PERTHESE™ Esthea Breast Implants EMP**

Vol.	Width	Height	Proj.	Catalog #
120 cc	9,2 cm	7,7 cm	3,6 cm	EMP120
145 cc	10,2 cm	8,5 cm	3,9 cm	EMP145
165 cc	10,7 cm	9,1 cm	4,0 cm	EMP165
200 cc	11,2 cm	9,6 cm	4,3 cm	EMP200
230 cc	11,7 cm	10,1 cm	4,5 cm	EMP230
265 cc	12,2 cm	10,6 cm	4,6 cm	EMP265
305 cc	12,7 cm	11,1 cm	4,9 cm	EMP305
345 cc	13,1 cm	11,6 cm	5,0 cm	EMP345
390 cc	13,6 cm	12,0 cm	5,2 cm	EMP390
435 cc	14,0 cm	12,5 cm	5,5 cm	EMP435
475 cc	14,6 cm	13,0 cm	5,7 cm	EMP475
535 cc	15,1 cm	13,3 cm	5,9 cm	EMP535
585 cc	15,5 cm	13,6 cm	6,1 cm	EMP585

**PERTHESE™ Esthea Breast Implants EFP**

Vol.	Width	Height	Proj.	Catalog #
140 cc	9,2 cm	8,5 cm	3,8 cm	EFP140
180 cc	10,2 cm	9,2 cm	4,1 cm	EFP180
210 cc	10,7 cm	9,7 cm	4,3 cm	EFP210
250 cc	11,2 cm	10,2 cm	4,5 cm	EFP250
285 cc	11,7 cm	10,7 cm	4,7 cm	EFP285
325 cc	12,2 cm	11,3 cm	4,9 cm	EFP325
370 cc	12,7 cm	11,8 cm	5,1 cm	EFP370
415 cc	13,1 cm	12,2 cm	5,3 cm	EFP415
470 cc	13,6 cm	12,6 cm	5,5 cm	EFP470
520 cc	14,0 cm	13,1 cm	5,7 cm	EFP520
575 cc	14,6 cm	13,5 cm	5,9 cm	EFP575
640 cc	15,1 cm	14,0 cm	6,1 cm	EFP640
710 cc	15,5 cm	14,2 cm	6,3 cm	EFP710

**PERTHESE™ Esthea Breast Implants EMHP**

Vol.	Width	Height	Proj.	Catalog #
195 cc	10,7 cm	9,1 cm	4,5 cm	EMHP195
235 cc	11,2 cm	9,6 cm	4,7 cm	EMHP235
270 cc	11,7 cm	10,1 cm	4,9 cm	EMHP270
320 cc	12,2 cm	10,6 cm	5,1 cm	EMHP320
360 cc	12,7 cm	11,1 cm	5,3 cm	EMHP360
405 cc	13,1 cm	11,6 cm	5,5 cm	EMHP405
455 cc	13,6 cm	12,0 cm	5,7 cm	EMHP455
490 cc	14,0 cm	12,5 cm	5,8 cm	EMHP490
545 cc	14,6 cm	13,0 cm	6,0 cm	EMHP545
625 cc	15,1 cm	13,3 cm	6,2 cm	EMHP625
680 cc	15,5 cm	13,6 cm	6,4 cm	EMHP680

**PERTHESE™ Esthea Breast Implants EMHP**


Vol.	Width	Height	Proj.	Catalog #
130 cc	9,2 cm	7,7 cm	4,0 cm	EMHP130
170 cc	10,2 cm	8,5 cm	4,3 cm	EMHP170
195 cc	10,7 cm	9,1 cm	4,5 cm	EMHP195
235 cc	11,2 cm	9,6 cm	4,7 cm	EMHP235
270 cc	11,7 cm	10,1 cm	4,9 cm	EMHP270
320 cc	12,2 cm	10,6 cm	5,1 cm	EMHP320
360 cc	12,7 cm	11,1 cm	5,3 cm	EMHP360
405 cc	13,1 cm	11,6 cm	5,5 cm	EMHP405
455 cc	13,6 cm	12,0 cm	5,7 cm	EMHP455
490 cc	14,0 cm	12,5 cm	5,8 cm	EMHP490
545 cc	14,6 cm	13,0 cm	6,0 cm	EMHP545
625 cc	15,1 cm	13,3 cm	6,2 cm	EMHP625
680 cc	15,5 cm	13,6 cm	6,4 cm	EMHP680

**PERTHESE™ Esthea Breast Implants AX**

Vol.	Width	Height	Proj.	Catalog #
190 cc	10,6 cm	10,8 cm	4,3 cm	AX190
210 cc	10,9 cm	11,2 cm	4,4 cm	AX210
240 cc	11,5 cm	12,2 cm	4,6 cm	AX240
260 cc	11,9 cm	12,5 cm	4,8 cm	AX260
285 cc	12,1 cm	12,8 cm	4,9 cm	AX285
300 cc	12,6 cm	13,2 cm	4,9 cm	AX300
340 cc	12,9 cm	13,7 cm	4,9 cm	AX340
380 cc	13,3 cm	14,1 cm	5,0 cm	AX380
400 cc	13,6 cm	14,2 cm	5,1 cm	AX400
450 cc	13,9 cm	14,7 cm	5,2 cm	AX450
480 cc	14,5 cm	15,2 cm	5,3 cm	AX480
500 cc	14,9 cm	15,6 cm	5,3 cm	AX500

**STAGE 2**

Selection of PERTHESE™ Esthea Implant to use as a long-term implant



**PERTHESE™ Esthea Breast Implants EFP**

Vol.	Width	Height	Proj.	Catalog #
140 cc	9,2 cm	8,5 cm	3,8 cm	EFP140
180 cc	10,2 cm	9,2 cm	4,1 cm	EFP180
210 cc	10,7 cm	9,7 cm	4,3 cm	EFP210
250 cc	11,2 cm	10,2 cm	4,5 cm	EFP250
285 cc	11,7 cm	10,7 cm	4,7 cm	EFP285
325 cc	12,2 cm	11,3 cm	4,9 cm	EFP325
370 cc	12,7 cm	11,8 cm	5,1 cm	EFP370
415 cc	13,1 cm	12,2 cm	5,3 cm	EFP415
470 cc	13,6 cm	12,6 cm	5,5 cm	EFP470
520 cc	14,0 cm	13,1 cm	5,7 cm	EFP520
575 cc	14,6 cm	13,5 cm	5,9 cm	EFP575
640 cc	15,1 cm	14,0 cm	6,1 cm	EFP640
710 cc	15,5 cm	14,2 cm	6,3 cm	EFP710

\*PLEASE NOTE: The MENTOR® CPX2™ Tissue Expander Injection Dome is not MRI compatible. For detailed indications, contraindications, warning and precautions associated with the use of this device, please refer to the Product Insert Data Sheet.

Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted.

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Silicone Gel Breast Implant Sizers, Round Shapes

#### PERTHESE™ Breast Implants TX, Sterile Gel Sizers\* Low Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
TX 100	100 cc	9,1 cm	2,2 cm	GABTX100STE
TX 125	125 cc	9,7 cm	2,4 cm	GABTX125STE
TX 150	150 cc	10,2 cm	2,6 cm	GABTX150STE
TX 170	170 cc	10,8 cm	2,7 cm	GABTX170STE
TX 185	185 cc	11,1 cm	2,8 cm	GABTX185STE
TX 200	200 cc	11,4 cm	2,9 cm	GABTX200STE
TX 220	220 cc	11,8 cm	2,9 cm	GABTX220STE
TX 240	240 cc	12,2 cm	3,0 cm	GABTX240STE
TX 260	260 cc	12,3 cm	3,1 cm	GABTX260STE
TX 280	280 cc	12,5 cm	3,2 cm	GABTX280STE
TX 300	300 cc	12,8 cm	3,2 cm	GABTX300STE
TX 325	325 cc	13,1 cm	3,3 cm	GABTX325STE
TX 350	350 cc	13,5 cm	3,4 cm	GABTX350STE

#### PERTHESE™ Breast Implants MX, Sterile Gel Sizers\* Moderate Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
MX 100	100 cc	8,7 cm	2,8 cm	GABMX100STE
MX 125	125 cc	9,2 cm	3,0 cm	GABMX125STE
MX 150	150 cc	9,9 cm	3,2 cm	GABMX150STE
MX 175	175 cc	10,4 cm	3,3 cm	GABMX175STE
MX 200	200 cc	11,0 cm	3,4 cm	GABMX200STE
MX 225	225 cc	11,5 cm	3,5 cm	GABMX225STE
MX 250	250 cc	12,0 cm	3,6 cm	GABMX250STE
MX 275	275 cc	12,2 cm	3,7 cm	GABMX275STE
MX 300	300 cc	12,5 cm	3,8 cm	GABMX300STE
MX 325	325 cc	12,8 cm	4,0 cm	GABMX325STE
MX 350	350 cc	13,0 cm	4,2 cm	GABMX350STE
MX 375	375 cc	13,2 cm	4,4 cm	GABMX375STE
MX 400	400 cc	13,5 cm	4,5 cm	GABMX400STE
MX 450	450 cc	13,9 cm	4,6 cm	GABMX450STE
MX 500	500 cc	14,2 cm	5,0 cm	GABMX500STE
MX 550	550 cc	14,5 cm	5,2 cm	GABMX550STE
MX 600	600 cc	14,9 cm	5,4 cm	GABMX600STE
MX 650	650 cc	15,2 cm	5,5 cm	GABMX650STE
MX 700	700 cc	15,6 cm	5,6 cm	GABMX700STE
MX 800	800 cc	16,1 cm	6,0 cm	GABMX800STE

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Silicone Gel Breast Implant Sizers, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 540, Sterile Gel Sizers\* High Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
540-075	75 cc	8,1 cm	2,3 cm	GAB540075STE
540-100	100 cc	8,2 cm	3,0 cm	GAB540100STE
540-125	125 cc	8,8 cm	3,2 cm	GAB540125STE
540-150	150 cc	9,2 cm	3,4 cm	GAB540150STE
540-175	175 cc	9,9 cm	3,6 cm	GAB540175STE
540-200	200 cc	10,1 cm	3,8 cm	GAB540200STE
540-230	230 cc	10,4 cm	4,0 cm	GAB540230STE
540-260	260 cc	11,0 cm	4,3 cm	GAB540260STE
540-280	280 cc	11,1 cm	4,4 cm	GAB540280STE
540-300	300 cc	11,2 cm	4,5 cm	GAB540300STE
540-325	325 cc	11,7 cm	4,6 cm	GAB540325STE
540-350	350 cc	11,9 cm	4,8 cm	GAB540350STE
540-380	380 cc	12,0 cm	4,9 cm	GAB540380STE
540-400	400 cc	12,1 cm	4,9 cm	GAB540400STE
540-440	440 cc	12,7 cm	5,3 cm	GAB540440STE
540-500	500 cc	13,2 cm	5,4 cm	GAB540500STE
540-550	550 cc	13,25 cm	5,9 cm	GAB540550STE
540-600	600 cc	13,8 cm	6,0 cm	GAB540600STE
540-650	650 cc	14,2 cm	6,1 cm	GAB540650STE
540-850	850 cc	15,5 cm	6,7 cm	GAB540850STE

\*Sterile - Single Use Only

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Inflatable Sizers, Round Shapes

#### PERTHESE™ Breast Implants TX, Sterile Inflatable Sizers\* Low Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
TX 100	100 cc	9,1 cm	2,2 cm	GGSTX100
TX 125	125 cc	9,7 cm	2,4 cm	GGSTX125
TX 150	150 cc	10,2 cm	2,6 cm	GGSTX150
TX 170	170 cc	10,8 cm	2,7 cm	GGSTX170
TX 185	185 cc	11,1 cm	2,8 cm	GGSTX185
TX 200	200 cc	11,4 cm	2,9 cm	GGSTX200
TX 220	220 cc	11,8 cm	2,9 cm	GGSTX220
TX 240	240 cc	12,2 cm	3,0 cm	GGSTX240
TX 260	260 cc	12,3 cm	3,1 cm	GGSTX260
TX 280	280 cc	12,5 cm	3,2 cm	GGSTX280
TX 300	300 cc	12,8 cm	3,2 cm	GGSTX300
TX 325	325 cc	13,1 cm	3,3 cm	GGSTX325
TX 350	350 cc	13,5 cm	3,4 cm	GGSTX350

#### PERTHESE™ Breast Implants MX, Sterile Inflatable Sizers\* Moderate Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
MX 125	125 cc	9,2 cm	3,0 cm	GGSMX125
MX 150	150 cc	9,9 cm	3,2 cm	GGSMX150
MX 175	175 cc	10,4 cm	3,3 cm	GGSMX175
MX 200	200 cc	11,0 cm	3,4 cm	GGSMX200
MX 225	225 cc	11,5 cm	3,5 cm	GGSMX225
MX 250	250 cc	12,0 cm	3,6 cm	GGSMX250
MX 275	275 cc	12,2 cm	3,7 cm	GGSMX275
MX 300	300 cc	12,5 cm	3,8 cm	GGSMX300
MX 325	325 cc	12,8 cm	4,0 cm	GGSMX325
MX 350	350 cc	13,0 cm	4,2 cm	GGSMX350
MX 375	375 cc	13,2 cm	4,4 cm	GGSMX375
MX 400	400 cc	13,5 cm	4,5 cm	GGSMX400

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Inflatable Sizers, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 540, Sterile Inflatable Sizers\* High Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
540-200	200 cc	10,1 cm	3,8 cm	GG540200
540-230	230 cc	10,4 cm	4,0 cm	GG540230
540-260	260 cc	11,0 cm	4,3 cm	GG540260
540-280	280 cc	11,1 cm	4,4 cm	GG540280
540-300	300 cc	11,2 cm	4,5 cm	GG540300
540-325	325 cc	11,7 cm	4,6 cm	GG540325
540-350	350 cc	11,9 cm	4,8 cm	GG540350
540-380	380 cc	12,0 cm	4,9 cm	GG540380
540-400	400 cc	12,1 cm	4,9 cm	GG540400

\*Sterile - Single Use Only

Sterile Silicone Gel Breast Implant Sizers,  
Anatomical Shapes

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*

ELP - EHP - ESHP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
ELP 190	190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	GABELP190STE
EHP 225	225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	GABEHP225STE
ESHP 235	235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	GABESHP235STE
ELP 240	240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	GABELP240STE
EHP 275	275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	GABEHP275STE
ESHP 295	295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	GABESHP295STE
ELP 270	270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	GABELP270STE
EHP 325	325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	GABEHP325STE
ESHP 355	355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	GABESHP355STE
ELP 310	310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	GABELP310STE
EHP 375	375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	GABEHP375STE
ESHP 415	415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	GABESHP415STE
ELP 360	360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	GABELP360STE
EHP 425	425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	GABEHP425STE
ESHP 475	475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	GABESHP475STE

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*  
AX, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
AX 190	190 cc	10,6 cm	10,8 cm	4,3 cm	3,7 cm	5,8 cm	7,8 cm	GABAX190STE
AX 210	210 cc	10,9 cm	11,2 cm	4,4 cm	3,8 cm	6,0 cm	8,2 cm	GABAX210STE
AX 240	240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	GABAX240STE
AX 260	260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	GABAX260STE
AX 285	285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	GABAX285STE
AX 300	300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	GABAX300STE
AX 340	340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	GABAX340STE
AX 380	380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	GABAX380STE
AX 400	400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	GABAX400STE
AX 450	450 cc	13,9 cm	14,7 cm	5,2 cm	4,8 cm	7,5 cm	10,8 cm	GABAX450STE
AX 480	480 cc	14,5 cm	15,2 cm	5,3 cm	4,9 cm	7,7 cm	11,2 cm	GABAX480STE
AX 500	500 cc	14,9 cm	15,6 cm	5,3 cm	5,0 cm	7,9 cm	11,5 cm	GABAX500STE

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Sterile Silicone Gel Breast Implant Sizers,  
Anatomical Shapes CONT'D

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*

EMP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
EMP 080	80 cc	8,1 cm	6,8 cm	3,2 cm	2,4 cm	3,7 cm	4,9 cm	GABEMP080STE
EMP 100	100 cc	8,6 cm	7,3 cm	3,5 cm	2,7 cm	4,2 cm	5,3 cm	GABEMP100STE
EMP 120	120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	GABEMP120STE
EMP 145	145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	GABEMP145STE
EMP 165	165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	GABEMP165STE
EMP 200	200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	GABEMP200STE
EMP 230	230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	GABEMP230STE
EMP 265	265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	GABEMP265STE
EMP 305	305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	GABEMP305STE
EMP 345	345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	GABEMP345STE
EMP 390	390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	GABEMP390STE
EMP 435	435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	GABEMP435STE
EMP 475	475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	GABEMP475STE
EMP 535	535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	GABEMP535STE
EMP 585	585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	GABEMP585STE

\*Sterile - Single Use Only



Sterile Inflatable Sizers,  
Anatomical Shapes

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*

ELP - EHP - ESHP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
ELP 190	190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	GGSELP190
EHP 225	225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	GGSEHP225
ESHP 235	235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	GGSESH235
ELP 240	240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	GGSELP240
EHP 275	275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	GGSEHP275
ESHP 295	295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	GGSESH295
ELP 270	270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	GGSELP270
EHP 325	325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	GGSEHP325
ESHP 355	355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	GGSESH355
ELP 310	310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	GGSELP310
EHP 375	375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	GGSEHP375
ESHP 415	415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	GGSESH415
ELP 360	360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	GGSELP360
EHP 425	425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	GGSEHP425
ESHP 475	475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	GGSESH475

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*

AX, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
AX 240	240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	GGSA240
AX 260	260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	GGSA260
AX 285	285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	GGSA285
AX 300	300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	GGSA300
AX 340	340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	GGSA340
AX 380	380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	GGSA380
AX 400	400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	GGSA400

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Sterile Inflatable Sizers,  
Anatomical Shapes CONT'D

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*

EMP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
EMP 120	120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	GGSEMP120
EMP 145	145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	GGSEMP145
EMP 165	165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	GGSEMP165
EMP 200	200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	GGSEMP200
EMP 230	230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	GGSEMP230
EMP 265	265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	GGSEMP265
EMP 305	305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	GGSEMP305
EMP 345	345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	GGSEMP345
EMP 390	390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	GGSEMP390
EMP 435	435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	GGSEMP435
EMP 475	475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	GGSEMP475
EMP 535	535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	GGSEMP535
EMP 585	585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	GGSEMP585

\*Sterile - Single Use Only

# PERTHESE™ SILICONE SHEETINGS

## Sterile Silicone Sheetings

### Dimensions 40 x 40 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	1	LP 300-1
Non reinforced	0.25 mm	1	LP 300-3

### Dimensions 50 x 50 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	10	LP 600-1
Non reinforced	0.25 mm	10	LP 600-3
Non reinforced	0.50 mm	10	LP 600-5

### Dimensions 60 x 100 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	1	LP 400-1
Non reinforced	0.25 mm	1	LP 400-3
Non reinforced	0.50 mm	1	LP 400-5
Non reinforced	1.00 mm	1	LP 400-7
Non reinforced	1.50 mm	1	LP 400-9

### Dimensions 70 x 100 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	5	LP 700-1
Non reinforced	0.25 mm	5	LP 700-3
Non reinforced	0.50 mm	5	LP 700-5
Non reinforced	1.00 mm	5	LP 700-7
Non reinforced	1.50 mm	5	LP 700-9
Non reinforced - extra firm	1.00 mm	5	LP 702-1

### Dimensions 90 x 100 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	1	LP 900-1
Non reinforced	0.25 mm	1	LP 900-3
Non reinforced	0.50 mm	1	LP 900-5
Non reinforced	1.00 mm	1	LP 900-7
Non reinforced	1.50 mm	1	LP 900-9

### Dimensions 150 x 200 mm

Description	Thickness	Packaging (units)	Catalog #
Non reinforced	0.125 mm	2	LP 800-1
Non reinforced	0.25 mm	2	LP 800-3
Non reinforced	0.50 mm	2	LP 800-5
Non reinforced	1.00 mm	2	LP 800-7
Non reinforced	1.50 mm	2	LP 800-9
Non reinforced - extra firm	1.00 mm	2	LP 802-1
Non reinforced - extra firm	2.00 mm	2	LP 802-3
Non reinforced - extra firm	3.00 mm	2	LP 802-5

### Corporate Headquarters

Santa Barbara, CA 93111 USA  
 Tel: +1 805 879 6000  
 Fax: +1 805 967 7108  
 www.mentorwllc.com

### Operations Europe, Middle East, Africa

#### Germany, Austria & Switzerland

Mentor Deutschland GmbH  
 Tel: +49 811 600 50 0  
 Fax: +49 811 600 50 11

#### Spain & Portugal

Mentor Medical Systems Iberica, S.L.  
 Tel: +34 91 562 2700  
 Fax: +34 91 728 3367

#### Italy, Greece & Israel

Johnson & Johnson Medical S.p.a.  
 Tel: +39 06 911 94 500  
 Fax: +39 02 644 10 563

#### United Kingdom & Ireland

Johnson & Johnson Medical Ltd  
 Tel: +44 (0)800 973 082  
 Fax: +44 (0)134 486 4122

#### France & Benelux

ETHICON France, S.A.S.  
 Tel: +33 1 55 00 22 45  
 Fax: +33 1 55 00 28 11

#### Russia & CIS

Johnson & Johnson LLC  
 Tel: +7 495 580 7777  
 Fax: +7 495 580 7878

#### Croatia, Macedonia, Serbia, Slovenia & Bulgaria

Johnson & Johnson S.E. d.o.o.  
 Tel: +386 1 40 11 800  
 Fax: +386 1 40 11 801

#### Czech Republic, Slovak Republic, Hungary & Romania

Johnson & Johnson s.r.o.  
 Tel: +420 227 012 371  
 Fax: +420 227 012 300

#### Poland, Lithuania, Latvia & Estonia

Johnson & Johnson Poland sp. z o.o.  
 Tel: +48 22 237 83 05  
 Fax: +48 22 237 83 53

#### Sweden, Denmark, Norway, Finland & Iceland

Johnson & Johnson AB  
 Tel: +46 8 626 2200  
 Fax: +46 8 754 5850

#### South Africa

Johnson & Johnson Medical (Pty) Ltd  
 Tel: +27 11 265 1135  
 Fax: +27 11 265 1114

#### Turkey

Johnson & Johnson Medical Turkey  
 Tel: +90 216 538 21 21  
 Fax: +90 216 538 21 30

#### Middle East, Egypt & Pakistan

Johnson & Johnson Middle East Inc  
 Tel: +971 4 429 72 00  
 Fax: +971 4 299 43 30

#### Maghreb & Sub Saharan Africa

Johnson & Johnson Professional Export  
 Tel: +44 (0)1506 594 500  
 Fax: +44 (0)1506 594 731

